




W SERIES CONTAINER SYSTEMS



INFOGRAPHICS



Lid: Perforated
Base: Closed



Lid: Safety Lid
Base: Perforated



Lid: Perforated
Base: Perforated



Lid: Safety Lid
Base: Closed



Lid: Closed
Base: Closed

There are 5 different types for W series containers



BAHADIR

explore / bahadir.com.tr

Length x Width x Height (mm)

Outside Dimensions
Inside Dimensions

585 x 280 x 105 mm
550 x 260 x 85 mm

Y 110.10W
Y 120.10W
Y 130.10W
Y 140.10W
Y 150.10W
Y 160.10W

585 x 280 x 135 mm
550 x 260 x 115 mm

Y 110.13W
Y 120.13W
Y 130.13W
Y 140.13W
Y 150.13W
Y 160.13W

585 x 280 x 155 mm
550 x 260 x 135 mm

Y 110.15W
Y 120.15W
Y 130.15W
Y 140.15W
Y 150.15W
Y 160.15W

585 x 280 x 205 mm
550 x 260 x 185 mm

Y 110.20W
Y 120.20W
Y 130.20W
Y 140.20W
Y 150.20W
Y 160.20W

585 x 280 x 265 mm
550 x 260 x 245 mm

Y 110.26W
Y 120.26W
Y 130.26W
Y 140.26W
Y 150.26W
Y 160.26W

Y 112.10W
Y 122.10W
Y 132.10W
Y 142.10W
Y 152.10W
Y 162.10W

Y 112.13W
Y 122.13W
Y 132.13W
Y 142.13W
Y 152.13W
Y 162.13W

Y 112.15W
Y 122.15W
Y 132.15W
Y 142.15W
Y 152.15W
Y 162.15W

Y 112.20W
Y 122.20W
Y 132.20W
Y 142.20W
Y 152.20W
Y 162.20W

Y 112.26W
Y 122.26W
Y 132.26W
Y 142.26W
Y 152.26W
Y 162.26W

585 x 280 x 105 mm
550 x 260 x 75 mm

Y 111.10W
Y 121.10W
Y 131.10W
Y 141.10W
Y 151.10W
Y 161.10W

585 x 280 x 135 mm
550 x 260 x 105 mm

Y 111.13W
Y 121.13W
Y 131.13W
Y 141.13W
Y 151.13W
Y 161.13W

585 x 280 x 155 mm
550 x 260 x 125 mm

Y 111.15W
Y 121.15W
Y 131.15W
Y 141.15W
Y 151.15W
Y 161.15W

585 x 280 x 205 mm
550 x 260 x 175 mm

Y 111.20W
Y 121.20W
Y 131.20W
Y 141.20W
Y 151.20W
Y 161.20W

585 x 280 x 265 mm
550 x 260 x 235 mm

Y 111.26W
Y 121.26W
Y 131.26W
Y 141.26W
Y 151.26W
Y 161.26W

Y 113.10W
Y 123.10W
Y 133.10W
Y 143.10W
Y 153.10W
Y 163.10W

Y 113.13W
Y 123.13W
Y 133.13W
Y 143.13W
Y 153.13W
Y 163.13W

Y 113.15W
Y 123.15W
Y 133.15W
Y 143.15W
Y 153.15W
Y 163.15W

Y 113.20W
Y 123.20W
Y 133.20W
Y 143.20W
Y 153.20W
Y 163.20W

Y 113.26W
Y 123.26W
Y 133.26W
Y 143.26W
Y 153.26W
Y 163.26W

Types & Colors



Lid Perforated
Lid Closed



Lid Safety Lid
Lid Closed



Lid Perforated
Body Perforated



Lid Safety Lid
Body Perforated

W Series Container Systems

Full Size (1/1 Size)

- They are made from anodized aluminum material.
- Container lids are perforated; the base can be both perforated or closed.
 - Container lids have additional safety lids.
- Lids are available in 6 different color options.

Base (Closed)



Outside Dimensions - Inside Dimensions

585 x 280 x 90 mm - 550 x 250 x 35 mm
585 x 280 x 120 mm - 550 x 260 x 115 mm
585 x 280 x 140 mm - 550 x 250 x 135 mm
585 x 280 x 190 mm - 550 x 260 x 185 mm
585 x 280 x 250 mm - 550 x 260 x 245 mm

Codes

Y 105.10W
Y 105.13W
Y 105.15W
Y 105.20W
Y 105.26W

Base (Perforated)



585 x 280 x 90 mm - 550 x 250 x 75 mm
585 x 280 x 120 mm - 550 x 260 x 105 mm
585 x 280 x 140 mm - 550 x 260 x 125 mm
585 x 280 x 190 mm - 550 x 260 x 175 mm
585 x 280 x 250 mm - 550 x 260 x 235 mm

Y 106.10W
Y 106.13W
Y 106.15W
Y 106.20W
Y 106.26W

Lid (Perforated)

585 x 280 mm

Codes & Colors



Y 100.01
Y 100.02
Y 100.03
Y 100.04
Y 100.05
Y 100.06

Safety Lid

585 x 280 mm

Codes & Colors




Y 103.01B
Y 103.02B
Y 103.03B
Y 103.04B
Y 103.05B
Y 103.06B



BAHADIR

explore / bahadir.com.tr



Length x Width x Height (mm)

Outside Dimensions
Inside Dimensions

465 x 280 x 105 mm
430 x 260 x 85 mm

Y 210.10W
Y 220.10W
Y 230.10W
Y 240.10W
Y 250.10W
Y 260.10W

465 x 280 x 135 mm
430 x 260 x 115 mm

Y 210.13W
Y 220.13W
Y 230.13W
Y 240.13W
Y 250.13W
Y 260.13W

465 x 280 x 155 mm
430 x 260 x 135 mm

Y 210.15W
Y 220.15W
Y 230.15W
Y 240.15W
Y 250.15W
Y 260.15W

465 x 280 x 205 mm
430 x 260 x 185 mm

Y 210.20W
Y 220.20W
Y 230.20W
Y 240.20W
Y 250.20W
Y 260.20W

465 x 280 x 265 mm
430 x 260 x 245 mm

Y 210.26W
Y 220.26W
Y 230.26W
Y 240.26W
Y 250.26W
Y 260.26W

Y 212.10W
Y 222.10W
Y 232.10W
Y 242.10W
Y 252.10W
Y 262.10W

Y 212.13W
Y 222.13W
Y 232.13W
Y 242.13W
Y 252.13W
Y 262.13W

Y 212.15W
Y 222.15W
Y 232.15W
Y 242.15W
Y 252.15W
Y 262.15W

Y 212.20W
Y 222.20W
Y 232.20W
Y 242.20W
Y 252.20W
Y 262.20W

Y 212.26W
Y 222.26W
Y 232.26W
Y 242.26W
Y 252.26W
Y 262.26W

465 x 280 x 105 mm
430 x 260 x 75 mm

Y 211.10W
Y 221.10W
Y 231.10W
Y 241.10W
Y 251.10W
Y 261.10W

465 x 280 x 135 mm
430 x 260 x 105 mm

Y 211.13W
Y 221.13W
Y 231.13W
Y 241.13W
Y 251.13W
Y 261.13W

465 x 280 x 155 mm
430 x 260 x 125 mm

Y 211.15W
Y 221.15W
Y 231.15W
Y 241.15W
Y 251.15W
Y 261.15W

465 x 280 x 205 mm
430 x 260 x 175 mm

Y 211.20W
Y 221.20W
Y 231.20W
Y 241.20W
Y 251.20W
Y 261.20W

465 x 280 x 265 mm
430 x 260 x 235 mm

Y 211.26W
Y 221.26W
Y 231.26W
Y 241.26W
Y 251.26W
Y 261.26W

Y 213.10W
Y 223.10W
Y 233.10W
Y 243.10W
Y 253.10W
Y 263.10W

Y 213.13W
Y 223.13W
Y 233.13W
Y 243.13W
Y 253.13W
Y 263.13W

Y 213.15W
Y 223.15W
Y 233.15W
Y 243.15W
Y 253.15W
Y 263.15W

Y 213.20W
Y 223.20W
Y 233.20W
Y 243.20W
Y 253.20W
Y 263.20W

Y 213.26W
Y 223.26W
Y 233.26W
Y 243.26W
Y 253.26W
Y 263.26W

Types &
Colors



Lid: Perforated
Type: Classic



Lid: Safety Lid
Type: Classic



Lid: Perforated
Type: Perforated



Lid: Safety Lid
Type: Perforated

W Series Container Systems

Middle Size (3/4 Size)

- They are made from anodized aluminum material.
- Container lids are perforated; the base can be both perforated or closed.
- Container lids have additional safety lids.
- Lids are available in 6 different color options.

Base (Closed)



Outside Dimensions - Inside Dimensions

465 x 280 x 90 mm - 430 x 260 x 85 mm
465 x 280 x 120 mm - 430 x 260 x 115 mm
465 x 280 x 140 mm - 430 x 260 x 135 mm
465 x 280 x 190 mm - 430 x 260 x 185 mm
465 x 280 x 250 mm - 430 x 260 x 245 mm

Codes

Y 205.10W
Y 205.13W
Y 205.15W
Y 205.20W
Y 205.26W

Base (Perforated)



465 x 280 x 90 mm - 430 x 260 x 75 mm
465 x 280 x 120 mm - 430 x 260 x 105 mm
465 x 280 x 140 mm - 430 x 260 x 125 mm
465 x 280 x 190 mm - 430 x 260 x 175 mm
465 x 280 x 250 mm - 430 x 260 x 235 mm

Y 206.10W
Y 206.13W
Y 206.15W
Y 206.20W
Y 206.26W

Lid (Perforated)

465 x 280 mm

Codes & Colors



Y 200.01
Y 200.02
Y 200.03
Y 200.04
Y 200.05
Y 200.06

Safety Lid

465 x 280 mm

Codes & Colors



Y 203.01B
Y 203.02B
Y 203.03B
Y 203.04B
Y 203.05B
Y 203.06B



BAHADIR

explore / bahadir.com.tr

Length x Width x Height (mm)

Types & Colors



Lid Perforated
Base Closed



Lid Safety Lin
Base Closed



Lid Perforated
Base Perforated



Lid Safety Lin
Base Perforated

Outside Dimensions
Inside Dimensions

295 x 280 x 105 mm
260 x 260 x 85 mm

Y 310.10W
Y 320.10W
Y 330.10W
Y 340.10W
Y 350.10W
Y 360.10W

295 x 280 x 135 mm
260 x 260 x 115 mm

Y 310.13W
Y 320.13W
Y 330.13W
Y 340.13W
Y 350.13W
Y 360.13W

295 x 280 x 155 mm
260 x 260 x 135 mm

Y 310.15W
Y 320.15W
Y 330.15W
Y 340.15W
Y 350.15W
Y 360.15W

295 x 280 x 205 mm
260 x 260 x 185 mm

Y 310.20W
Y 320.20W
Y 330.20W
Y 340.20W
Y 350.20W
Y 360.20W

295 x 280 x 265 mm
260 x 260 x 245 mm

Y 310.26W
Y 320.26W
Y 330.26W
Y 340.26W
Y 350.26W
Y 360.26W

295 x 280 x 105 mm
260 x 260 x 75 mm

Y 311.10W
Y 321.10W
Y 331.10W
Y 341.10W
Y 351.10W
Y 361.10W

295 x 280 x 135 mm
260 x 260 x 105 mm

Y 311.13W
Y 321.13W
Y 331.13W
Y 341.13W
Y 351.13W
Y 361.13W

295 x 280 x 155 mm
260 x 260 x 125 mm

Y 311.15W
Y 321.15W
Y 331.15W
Y 341.15W
Y 351.15W
Y 361.15W

295 x 280 x 205 mm
260 x 260 x 175 mm

Y 311.20W
Y 321.20W
Y 331.20W
Y 341.20W
Y 351.20W
Y 361.20W

295 x 280 x 265 mm
260 x 260 x 235 mm

Y 311.26W
Y 321.26W
Y 331.26W
Y 341.26W
Y 351.26W
Y 361.26W

Y 313.10W
Y 323.10W
Y 333.10W
Y 343.10W
Y 353.10W
Y 363.10W

Y 313.13W
Y 323.13W
Y 333.13W
Y 343.13W
Y 353.13W
Y 363.13W

Y 313.15W
Y 323.15W
Y 333.15W
Y 343.15W
Y 353.15W
Y 363.15W

Y 313.20W
Y 323.20W
Y 333.20W
Y 343.20W
Y 353.20W
Y 363.20W

Y 313.26W
Y 323.26W
Y 333.26W
Y 343.26W
Y 353.26W
Y 363.26W

W Series Container Systems

Half Size (1/2 Size)

- They are made from anodized aluminum material.
- Container lids are perforated; the base can be both perforated or closed.
- Container lids have additional safety lids.
- Lids are available in 6 different color options.

Base (Closed)



Outside Dimensions - Inside Dimensions

295 x 280 x 90 mm - 260 x 260 x 85 mm
295 x 280 x 120 mm - 260 x 260 x 115 mm
295 x 280 x 140 mm - 260 x 260 x 135 mm
295 x 280 x 190 mm - 260 x 260 x 185 mm
295 x 280 x 250 mm - 260 x 260 x 245 mm

Codes

Y 305.10W
Y 305.13W
Y 305.15W
Y 305.20W
Y 305.26W

Base (Perforated)



295 x 280 x 90 mm - 260 x 260 x 75 mm
295 x 280 x 120 mm - 260 x 260 x 105 mm
295 x 280 x 140 mm - 260 x 260 x 125 mm
295 x 280 x 190 mm - 260 x 260 x 175 mm
295 x 280 x 250 mm - 260 x 260 x 235 mm

Y 306.10W
Y 306.13W
Y 306.15W
Y 306.20W
Y 306.26W

Lid (Perforated)

295 x 280 mm

Codes & Colors



Y 300.01
Y 300.02
Y 300.03
Y 300.04
Y 300.05
Y 300.06

Safety Lid

295 x 280 mm

Codes & Colors



Y 303.01B
Y 303.02B
Y 303.03B
Y 303.04B
Y 303.05B
Y 303.06B



W Series Container Systems

Wide Body

- They are made from anodized aluminum material.
- Container lids are perforated; the base can be both perforated or closed.
- Wide Body Containers provide ultimate efficiency with high sterilization capacity.
- Lids are available in 6 different color options.



BAHADIR

explore / bahadir.com.tr

Length x Width x Height (mm)




Outside Dimensions
Inside Dimensions

600 x 400 x 125 mm
565 x 380 x 105 mm

Y 110.62W
Y 120.62W
Y 130.62W
Y 140.62W
Y 150.62W
Y 160.62W

600 x 400 x 185 mm
565 x 380 x 165 mm

Y 110.68W
Y 120.68W
Y 130.68W
Y 140.68W
Y 150.68W
Y 160.68W



Outside Dimensions
Inside Dimensions

600 x 400 x 125 mm
565 x 380 x 95 mm

Y 111.62W
Y 121.62W
Y 131.62W
Y 141.62W
Y 151.62W
Y 161.62W

600 x 400 x 185 mm
565 x 380 x 155 mm

Y 111.68W
Y 121.68W
Y 131.68W
Y 141.68W
Y 151.68W
Y 161.68W

Base (Closed)



Outside Dimensions - Inside Dimensions

600 x 400 x 110 mm - 565 x 380 x 105 mm
600 x 400 x 170 mm - 565 x 380 x 165 mm

Codes

Y 105.62W
Y 105.68W

Lid (Perforated)

600 x 400 mm



Base (Perforated)



600 x 400 x 110 mm - 565 x 380 x 95 mm
600 x 400 x 170 mm - 565 x 380 x 155 mm

Y 106.62W
Y 106.68W

Codes & Colors

Y 100.51
Y 100.52
Y 100.53
Y 100.54
Y 100.55
Y 100.56





W Series Container Systems Long Body

- They are made from anodized aluminum material.
- Container lids are perforated; the base can be both perforated or closed.
- Long Body Containers provide ultimate efficiency with high sterilization capacity.
- Lids are available in 6 different color options.



Outside Dimensions
Inside Dimensions

710 x 280 x 205 mm
685 x 260 x 185 mm

Y 110.70W
Y 120.70W
Y 130.70W
Y 140.70W
Y 150.70W
Y 160.70W



Outside Dimensions
Inside Dimensions

710 x 280 x 205 mm
685 x 260 x 175 mm

Y 111.70W
Y 121.70W
Y 131.70W
Y 141.70W
Y 151.70W
Y 161.70W

Base (Closed)



Outside Dimensions - Inside Dimensions

710 x 280 x 190 mm - 685 x 260 x 185 mm

Codes

Y 105.70W

Lid (Perforated)

710 x 280 mm



Base (Perforated)



710 x 280 x 190 mm - 685 x 260 x 175 mm

Y 106.70W

Codes & Colors

Y 100.71
Y 100.72
Y 100.73
Y 100.74
Y 100.75
Y 100.76





W Series Container Systems

Removal Full Size (1/1 Size)

They are made from stainless steel material.
Removal Containers are used for the transport and disinfection
of medical waste.



explore / bahadir.com.tr

Length x Width x Height (mm)



Outside Dimensions
Inside Dimensions

585 x 280 x 155 mm
550 x 260 x 135 mm

Y 115.15W
Y 125.15W
Y 135.15W
Y 145.15W
Y 155.15W
Y 165.15W

Base (Closed)



Outside Dimensions - Inside Dimensions

585 x 280 x 140 mm - 550 x 260 x 135 mm

Codes

Y 105.15W

Lid (Closed)

585 x 280 mm



Codes & Colors

Y 102.01
Y 102.02
Y 102.03
Y 102.04
Y 102.05
Y 102.06





W Series Container Systems

Removal Half Size (1/2 Size)

- They are made from stainless steel material
- Removal Containers are used for the transport and disinfection of medical waste



Outside Dimensions Inside Dimensions

295 x 280 x 105 mm
260 x 260 x 85 mm

Y 315.10W
Y 325.10W
Y 335.10W
Y 345.10W
Y 355.10W
Y 365.10W

295 x 280 x 155 mm
260 x 260 x 135 mm

Y 315.15W
Y 325.15W
Y 335.15W
Y 345.15W
Y 355.15W
Y 365.15W

Base (Closed)



Outside Dimensions - Inside Dimensions

295 x 280 x 90 mm - 260 x 260 x 85 mm
295 x 280 x 140 mm - 260 x 260 x 135 mm

Codes

Y 305.10W
Y 305.15W

Lid (Closed)

295 x 280 mm

Codes & Colors



Y 302.01
Y 302.02
Y 302.03
Y 302.04
Y 302.05
Y 302.06





A SERIES CONTAINER SYSTEMS



INFOGRAPHICS



Lid: Perforated
Base: Closed



Lid: Perforated
Base: Perforated

There are 2 different types for A series containers



A Series Container Systems Flat Model

- They are made from anodized aluminum material.
- Lids are available in 6 different color options.



BAHADIR

explore / bahadir.com.tr

Length x Width x Height (mm)



Outside Dimensions
Inside Dimensions

300 x 280 x 55 mm
270 x 270 x 45 mm

Y 310.05A
Y 320.05A
Y 330.05A
Y 340.05A
Y 350.05A
Y 360.05A

300 x 280 x 85 mm
270 x 270 x 75 mm

Y 310.08A
Y 320.08A
Y 330.08A
Y 340.08A
Y 350.08A
Y 360.08A



Outside Dimensions
Inside Dimensions

300 x 280 x 85 mm
270 x 270 x 75 mm

Y 311.08A
Y 321.08A
Y 331.08A
Y 341.08A
Y 351.08A
Y 361.08A

Base (Closed)



Outside Dimensions - Inside Dimensions

300 x 280 x 45 mm - 270 x 270 x 45 mm
300 x 280 x 75 mm - 270 x 270 x 75 mm

Codes

Y 305.05A
Y 305.08A

Lid (Perforated)

300 x 280 mm



Base (Perforated)



300 x 280 x 75 mm - 270 x 270 x 65 mm

Y 306.08A

Codes & Colors

Y 304.01A
Y 304.02A
Y 304.03A
Y 304.04A
Y 304.05A
Y 304.06A



A Series Container Systems Mini Model

- They are made from anodized aluminum material.
- Type and sizes are suitable for the sterilization of micro instruments.
- Lids are available in 6 different color options.



Outside Dimensions
Inside Dimensions

310 x 140 x 45 mm
280 x 130 x 40 mm

Y 410.04A
Y 420.04A
Y 430.04A
Y 440.04A
Y 450.04A
Y 460.04A

310 x 140 x 75 mm
280 x 130 x 70 mm

Y 410.07A
Y 420.07A
Y 430.07A
Y 440.07A
Y 450.07A
Y 460.07A

310 x 140 x 105 mm
280 x 130 x 100 mm

Y 410.10A
Y 420.10A
Y 430.10A
Y 440.10A
Y 450.10A
Y 460.10A

310 x 140 x 45 mm
280 x 130 x 35 mm

Y 411.04A
Y 421.04A
Y 431.04A
Y 441.04A
Y 451.04A
Y 461.04A

310 x 140 x 75 mm
280 x 130 x 65 mm

Y 411.07A
Y 421.07A
Y 431.07A
Y 441.07A
Y 451.07A
Y 461.07A

310 x 140 x 105 mm
280 x 130 x 95 mm

Y 411.10A
Y 421.10A
Y 431.07A
Y 441.10A
Y 451.10A
Y 461.10A

Base (Closed)



Outside Dimensions - Inside Dimensions

310 x 130 x 40 mm - 280 x 130 x 40 mm
310 x 130 x 70 mm - 280 x 130 x 70 mm
310 x 130 x 100 mm - 280 x 130 x 100 mm

Codes

Y 405.04A
Y 405.07A
Y 405.10A

Lid (Perforated)



310 x 140 mm

Base (Perforated)



310 x 130 x 40 mm - 280 x 130 x 35 mm
310 x 130 x 70 mm - 280 x 130 x 65 mm

Y 406.04A
Y 406.07A

Codes & Colors

Y 400.01A
Y 400.02A
Y 400.03A
Y 400.04A



A Series Container Systems Dental Model

- They are made from anodized aluminum material.
- Type and sizes are suitable for the sterilization of dental instruments.
- Lids are available in 6 different color options.



BAHADIR

explore / bahadir.com.tr

Length x Width x Height (mm)

Outside Dimensions Inside Dimensions

310 x 190 x 45 mm 280 x 180 x 40 mm	310 x 190 x 65 mm 280 x 180 x 60 mm	310 x 190 x 85 mm 280 x 180 x 80 mm	310 x 190 x 105 mm 280 x 180 x 100 mm	310 x 190 x 135 mm 280 x 180 x 130 mm
Y 510.04A	Y 510.06A	Y 510.08A	Y 510.10A	Y 510.13A
Y 520.04A	Y 520.06A	Y 520.08A	Y 520.10A	Y 520.13A
Y 530.04A	Y 530.06A	Y 530.08A	Y 530.10A	Y 530.13A
Y 540.04A	Y 540.06A	Y 540.08A	Y 540.10A	Y 540.13A
Y 550.04A	Y 550.06A	Y 550.08A	Y 550.10A	Y 550.13A
Y 560.04A	Y 560.06A	Y 560.08A	Y 560.10A	Y 560.13A
310 x 190 x 45 mm 280 x 180 x 35 mm	310 x 190 x 65 mm 280 x 180 x 55 mm	310 x 190 x 85 mm 280 x 180 x 75 mm	310 x 190 x 105 mm 280 x 180 x 95 mm	310 x 190 x 135 mm 280 x 180 x 125 mm
Y 511.04A	Y 511.06A	Y 511.08A	Y 511.10A	Y 511.13A
Y 521.04A	Y 521.06A	Y 521.08A	Y 521.10A	Y 521.13A
Y 531.04A	Y 531.06A	Y 531.08A	Y 531.10A	Y 531.13A
Y 541.04A	Y 541.06A	Y 541.08A	Y 541.10A	Y 541.13A
Y 551.04A	Y 551.06A	Y 551.08A	Y 551.10A	Y 551.13A
Y 561.04A	Y 561.06A	Y 561.08A	Y 561.10A	Y 561.13A

Base (Closed)

Outside Dimensions - Inside Dimensions

Codes


Lid (Perforated)

310 x 190 mm



310 x 190 x 40 mm - 280 x 180 x 40 mm
310 x 190 x 60 mm - 280 x 180 x 60 mm
310 x 190 x 80 mm - 280 x 180 x 80 mm
310 x 190 x 100 mm - 280 x 180 x 100 mm
310 x 190 x 130 mm - 280 x 180 x 130 mm

Y 505.04A
Y 505.06A
Y 505.08A
Y 505.10A
Y 505.13A



Base (Perforated)



310 x 190 x 40 mm - 280 x 180 x 35 mm
310 x 190 x 60 mm - 280 x 180 x 55 mm
310 x 190 x 80 mm - 280 x 180 x 75 mm

Y 506.04A
Y 506.06A
Y 506.08A

Codes & Colors

Y 500.01A
Y 500.02A
Y 500.03A





A Series Container Systems

Mini Dental Model

- They are made from anodized aluminum material.
- Type and sizes are suitable for the sterilization of micro dental instruments.
- Lids are available in 6 different color options.



Outside Dimensions
Inside Dimensions

210 x 155 x 55 mm
180 x 145 x 50 mm

Y 515.05A
Y 525.05A
Y 535.05A
Y 545.05A
Y 555.05A
Y 565.05A

210 x 155 x 65 mm
180 x 145 x 60 mm

Y 515.06A
Y 525.06A
Y 535.06A
Y 545.06A
Y 555.06A
Y 565.06A

Base (Closed)



Outside Dimensions - Inside Dimensions

210 x 155 x 50 mm - 180 x 145 x 50 mm
210 x 155 x 60 mm - 180 x 145 x 60 mm

Codes

Y 508.05A
Y 508.06A

Lid (Perforated)

210 x 155 mm



Codes & Colors

Y 501.01A
Y 501.02A
Y 501.03A
Y 501.04A
Y 501.05A
Y 501.06A



E Series - Dental Model

68

E Series - Mini Dental Model

69

E Series - Flat Model

70

E Series - Mini Model

71

pages

CONTENTS



E SERIES CONTAINER SYSTEMS



INFOGRAPHICS



Lid: Perforated
Base: Closed



Lid: Perforated
Base: Perforated

There are 5 different types for E series containers



E Series Container Systems

Dental Model

- They are made from anodized aluminum material.
- Type and sizes are suitable for the sterilization of dental instruments.
- Lids are available in 6 different color options.



explore / bahadir.com.tr

Length x Width x Height (mm)



Outside Dimensions
Inside Dimensions

Outside Dimensions	Inside Dimensions	Outside Dimensions	Inside Dimensions	Outside Dimensions	Inside Dimensions	Outside Dimensions	Inside Dimensions	Outside Dimensions	Inside Dimensions
310 x 190 x 45 mm 280 x 180 x 40 mm	Y 510.04E Y 520.04E Y 530.04E Y 540.04E Y 550.04E Y 560.04E	310 x 190 x 65 mm 280 x 180 x 60 mm	Y 510.06E Y 520.06E Y 530.06E Y 540.06E Y 550.06E Y 560.06E	310 x 190 x 85 mm 280 x 180 x 80 mm	Y 510.08E Y 520.08E Y 530.08E Y 540.08E Y 550.08E Y 560.08E	310 x 190 x 105 mm 280 x 180 x 100 mm	Y 510.10E Y 520.10E Y 530.10E Y 540.10E Y 550.10E Y 560.10E	310 x 190 x 135 mm 280 x 180 x 130 mm	Y 510.13E Y 520.13E Y 530.13E Y 540.13E Y 550.13E Y 560.13E
310 x 190 x 45 mm 280 x 180 x 35 mm	Y 511.04E Y 521.04E Y 531.04E Y 541.04E Y 551.04E Y 561.04E	310 x 190 x 65 mm 280 x 180 x 55 mm	Y 511.06E Y 521.06E Y 531.06E Y 541.06E Y 551.06E Y 561.06E	310 x 190 x 85 mm 280 x 180 x 75 mm	Y 511.08E Y 521.08E Y 531.08E Y 541.08E Y 551.08E Y 561.08E	310 x 190 x 105 mm 280 x 180 x 95 mm	Y 511.10E Y 521.10E Y 531.10E Y 541.10E Y 551.10E Y 561.10E	310 x 190 x 135 mm 280 x 180 x 125 mm	Y 511.13E Y 521.13E Y 531.13E Y 541.13E Y 551.13E Y 561.13E

Base (Closed)

Outside Dimensions - Inside Dimensions

Codes

Lid (Perforated)

310 x 190 mm



310 x 190 x 40 mm - 280 x 180 x 40 mm
310 x 190 x 60 mm - 280 x 180 x 60 mm
310 x 190 x 80 mm - 280 x 180 x 80 mm
310 x 190 x 100 mm - 280 x 180 x 100 mm
310 x 190 x 130 mm - 280 x 180 x 130 mm

Y 505.04E
Y 505.06E
Y 505.08E
Y 505.10E
Y 505.13E



Base (Perforated)

310 x 190 x 40 mm - 280 x 180 x 35 mm
310 x 190 x 60 mm - 280 x 180 x 55 mm
310 x 190 x 80 mm - 280 x 180 x 75 mm

Y 506.04E
Y 506.06E
Y 506.08E

Codes & Colors

Y 500.01E
Y 500.02E
Y 500.03E





E Series Container Systems Mini Dental Model

- They are made from anodized aluminum material.
- Type and sizes are suitable for the sterilization of micro dental instruments.
- Lids are available in 6 different color options.



Outside Dimensions
Inside Dimensions

210 x 155 x 55 mm
180 x 145 x 50 mm

Y 515.05E
Y 525.05E
Y 535.05E
Y 545.05E
Y 555.05E
Y 565.05E

210 x 155 x 65 mm
180 x 145 x 60 mm

Y 515.06E
Y 525.06E
Y 535.06E
Y 545.06E
Y 555.06E
Y 565.06E



Outside Dimensions
Inside Dimensions

210 x 155 x 55 mm
180 x 145 x 50 mm

Y 516.05E
Y 526.05E
Y 536.05E
Y 546.05E
Y 556.05E
Y 566.05E

210 x 155 x 65 mm
180 x 145 x 60 mm

Y 516.06E
Y 526.06E
Y 536.06E
Y 546.06E
Y 556.06E
Y 566.06E

Base (Closed)

Outside Dimensions - Inside Dimensions

Codes

Lid (Perforated)

210 x 155 mm



210 x 155 x 50 mm - 180 x 145 x 50 mm
210 x 155 x 60 mm - 180 x 145 x 60 mm

Y 508.05E
Y 508.06E



Base (Perforated)

Outside Dimensions - Inside Dimensions

Codes

Codes & Colors



210 x 155 x 50 mm - 180 x 145 x 45 mm
210 x 155 x 60 mm - 180 x 145 x 55 mm

Y 509.05E
Y 509.06E

Y 501.01E
Y 501.02E
Y 501.03E
Y 501.04E
Y 501.05E
Y 501.06E



E Series Container Systems Flat Model

- They are made from anodized aluminum material.
- Lids are available in 6 different color options.



BAHADIR

explore / bahadir.com.tr

Length x Width x Height (mm)



Outside Dimensions
Inside Dimensions

300 x 280 x 55 mm
270 x 270 x 45 mm

Y 310.05E
Y 320.05E
Y 330.05E
Y 340.05E
Y 350.05E
Y 360.05E

300 x 280 x 85 mm
270 x 270 x 75 mm

Y 310.08E
Y 320.08E
Y 330.08E
Y 340.08E
Y 350.08E
Y 360.08E



Outside Dimensions
Inside Dimensions

300 x 280 x 85 mm
270 x 270 x 65 mm

Y 311.08E
Y 321.08E
Y 331.08E
Y 341.08E
Y 351.08E
Y 361.08E

Base (Closed)



Outside Dimensions - Inside Dimensions

300 x 280 x 45 mm - 270 x 270 x 45 mm
300 x 280 x 75 mm - 270 x 270 x 75 mm

Codes

Y 305.05E
Y 305.08E

Lid (Perforated)

300 x 280 mm



Base (Perforated)



300 x 280 x 75 mm - 270 x 270 x 65 mm

Y 306.08E

Codes & Colors

Y 304.01E
Y 304.02E
Y 304.03E
Y 304.04E
Y 304.05E
Y 304.06E



E Series Container Systems

Mini Model

- They are made from anodized aluminum material
- Type and sizes are suitable for the sterilization of dental instruments
- Lids are available in 6 different color options



Outside Dimensions
Inside Dimensions

310 x 140 x 45 mm
280 x 130 x 40 mm

Y 410.04E
Y 420.04E
Y 430.04E
Y 440.04E
Y 450.04E
Y 460.04E

310 x 140 x 75 mm
280 x 130 x 70 mm

Y 410.07E
Y 420.07E
Y 430.07E
Y 440.07E
Y 450.07E
Y 460.07E

310 x 140 x 105 mm
280 x 130 x 100 mm

Y 410.10E
Y 420.10E
Y 430.10E
Y 440.10E
Y 450.10E
Y 460.10E

310 x 140 x 45 mm
280 x 130 x 35 mm

Y 411.04E
Y 421.04E
Y 431.04E
Y 441.04E
Y 451.04E
Y 461.04E

310 x 140 x 75 mm
280 x 130 x 65 mm

Y 411.07E
Y 421.07E
Y 431.07E
Y 441.07E
Y 451.07E
Y 461.07E

310 x 140 x 105 mm
280 x 130 x 95 mm

Y 411.10E
Y 421.10E
Y 431.07E
Y 441.10E
Y 451.10E
Y 461.10E

Base (Closed)



Outside Dimensions - Inside Dimensions

310 x 130 x 40 mm - 280 x 130 x 40 mm
310 x 130 x 70 mm - 280 x 130 x 70 mm
310 x 130 x 100 mm - 280 x 130 x 100 mm

Codes

Y 405.04E
Y 405.07E
Y 405.10E

Lid (Perforated)

310 x 140 mm



Base (Perforated)



310 x 130 x 40 mm - 280 x 130 x 35 mm
310 x 130 x 70 mm - 280 x 130 x 65 mm

Y 406.04E
Y 406.07E

Codes & Colors

Y 400.01E
Y 400.02E
Y 400.03E



NEXT PAGE

PERFORATED TRAYS & WIRE BASKETS



Perforated Trays


Contents

	pages
Standard Model - Perforated Tray	82
Dental, Mini & Optic Model - Perforated Tray	83



BAHADIR

explore / bahadir.com.tr



Our design provides easy handling and carrying even the lid is closed.



With silicone holders and dividers, the organization of surgical instruments is easy and safe.



Mesh Dimensions: 4x4 mm
Range Of Mesh: 1 mm



Polymer foot optimizes stacking of perforated trays into the sterilization containers.



Polymer Foot - Y 835,00

Outside Dimensions ►

Length x Width x Height (mm)

Types



Basic



Basic Grid



Basic (With Feet)



Basic Grid (With Feet)

Half Size (1/2 Size)

Y 380.03
242 x 252 x 25 mm

Y 380.05
242 x 252 x 45 mm

Y 380.07
242 x 252 x 65 mm

Y 380.10
242 x 252 x 95 mm

Y 380.03L
246 x 256 x 25 mm

Y 380.05L
246 x 256 x 45 mm

Y 380.07L
246 x 256 x 65 mm

Y 380.10L
246 x 256 x 95 mm

Y 380.03F
242 x 252 x 40 mm

Y 380.05F
242 x 252 x 60 mm

Y 380.07F
242 x 252 x 80 mm

Y 380.10F
242 x 252 x 110 mm

Y 380.03FL
246 x 256 x 40 mm

Y 380.05FL
246 x 256 x 60 mm

Y 380.07FL
246 x 256 x 80 mm

Y 380.10FL
246 x 256 x 110 mm

Y 380.01

Middle Size (3/4 Size)

Y 280.03
405 x 252 x 25 mm

Y 280.05
405 x 252 x 45 mm

Y 280.07
405 x 252 x 65 mm

Y 280.10
405 x 252 x 95 mm

Y 280.03L
409 x 256 x 25 mm

Y 280.05L
409 x 256 x 45 mm

Y 280.07L
409 x 256 x 65 mm

Y 280.10L
409 x 256 x 95 mm

Y 280.03F
405 x 252 x 40 mm

Y 280.05F
405 x 252 x 60 mm

Y 280.07F
405 x 252 x 80 mm

Y 280.10F
405 x 252 x 110 mm

Y 280.03FL
409 x 256 x 40 mm

Y 280.05FL
409 x 256 x 60 mm

Y 280.07FL
409 x 256 x 80 mm

Y 280.10FL
409 x 256 x 110 mm

Y 280.01

DIN Size

Y 180.03
484 x 252 x 25 mm

Y 180.05
484 x 252 x 45 mm

Y 180.07
484 x 252 x 65 mm

Y 180.10
484 x 252 x 95 mm

Y 180.03L
488 x 256 x 25 mm

Y 180.05L
488 x 256 x 45 mm

Y 180.07L
488 x 256 x 65 mm

Y 180.10L
488 x 256 x 95 mm

Y 180.03F
484 x 252 x 40 mm

Y 180.05F
484 x 252 x 60 mm

Y 180.07F
484 x 252 x 80 mm

Y 180.10F
484 x 252 x 110 mm

Y 180.03FL
488 x 256 x 40 mm

Y 180.05FL
488 x 256 x 60 mm

Y 180.07FL
488 x 256 x 80 mm

Y 180.10FL
488 x 256 x 110 mm

Y 180.01

Full Size (1/1 Size)

Y 181.03
539 x 252 x 25 mm

Y 181.05
539 x 252 x 45 mm

Y 181.07
539 x 252 x 65 mm

Y 181.10
539 x 252 x 95 mm

Y 181.03L
543 x 256 x 25 mm

Y 181.05L
543 x 256 x 45 mm

Y 181.07L
543 x 256 x 65 mm

Y 181.10L
543 x 256 x 95 mm

Y 181.03F
539 x 252 x 40 mm

Y 181.05F
539 x 252 x 60 mm

Y 181.07F
539 x 252 x 80 mm

Y 181.10F
539 x 252 x 110 mm

Y 181.03FL
543 x 256 x 40 mm

Y 181.05FL
543 x 256 x 60 mm

Y 181.07FL
543 x 256 x 80 mm


Y 181.10FL
543 x 256 x 110 mm

Y 181.01

Perforated Trays Optic Model



Outside Dimensions



Single Scope

Y 695.00
272 x 83 x 44 mm

Y 696.00
452 x 83 x 44 mm

Outside Dimensions



Double Scope

Y 695.03
272 x 103 x 44 mm

Y 696.03
452 x 103 x 44 mm

Perforated Trays Mini Model & Dental Model



Outside Dimensions



Mini Model

Y 480.04
270 x 123 x 45 mm

Y 480.06
270 x 123 x 65 mm

Outside Dimensions



Dental Model

Y 580.05
270 x 170 x 50 mm

Y 580.10
270 x 170 x 100 mm

Hybrid Baskets


Contents

	pages
Standard Hybrid - Intermediate Crimp	86
Standard Hybrid - Plain Weave	87
Standard Hybrid - Welded Wire Mesh	88
Specific Hybrid - Plain Weave	89



BAHADIR

explore / bahadir.com.tr



Base: Intermediate Crimp
Mesh Dimensions: 5x5 mm
Wire Diameter: 1 mm
Sides: Perforated sheet
Perforation Diameter: 3 mm



Base: Plain Weave
Mesh Dimensions: 4x4 mm
Wire Diameter: 1 mm
Sides: Perforated sheet
Perforation Diameter: 3 mm




Base: Welded Wire Mesh
Mesh Dimensions: 6x6 mm
Wire Diameter: 1,2 mm
Sides: Perforated sheet
Perforation Diameter: 3 mm



BAHADIR

explore / bahadir.com.tr

Standard Hybrid Baskets Intermediate Crimp




Length x Width x Height (mm)

Outside Dimensions ▶

Types ▼ Page Date to End 13	Half Size (1/2 Size)	Middle Size (3/4 Size)	DIN Size	Full Size (1/1 Size)
	Y 396.05A 245 x 255 x 50 mm	Y 296.05A 405 x 255 x 50 mm	Y 196.05A 485 x 255 x 50 mm	Y 197.05A 540 x 255 x 50 mm
	Y 396.06A 245 x 255 x 60 mm	Y 296.06A 405 x 255 x 60 mm	Y 196.06A 485 x 255 x 60 mm	Y 197.06A 540 x 255 x 60 mm
	Y 396.07A 245 x 255 x 70 mm	Y 296.07A 405 x 255 x 70 mm	Y 196.07A 485 x 255 x 70 mm	Y 197.07A 540 x 255 x 70 mm
	Y 396.10A 245 x 255 x 100 mm	Y 296.10A 405 x 255 x 100 mm	Y 196.10A 485 x 255 x 100 mm	Y 197.10A 540 x 255 x 100 mm
	Y 396.05AL 245 x 255 x 50 mm	Y 296.05AL 405 x 255 x 50 mm	Y 196.05AL 485 x 255 x 50 mm	Y 197.05AL 540 x 255 x 50 mm
	Y 396.06AL 245 x 255 x 60 mm	Y 296.06AL 405 x 255 x 60 mm	Y 196.06AL 485 x 255 x 60 mm	Y 197.06AL 540 x 255 x 60 mm
	Y 396.07AL 245 x 255 x 70 mm	Y 296.07AL 405 x 255 x 70 mm	Y 196.07AL 485 x 255 x 70 mm	Y 197.07AL 540 x 255 x 70 mm
	Y 396.10AL 245 x 255 x 100 mm	Y 296.10AL 405 x 255 x 100 mm	Y 196.10AL 485 x 255 x 100 mm	Y 197.10AL 540 x 255 x 100 mm
	Y 390.01 215 x 240 mm	Y 290.01 375 x 240 mm	Y 190.01 455 x 240 mm	Y 191.01 510 x 240 mm


Product Details

Base & Lid
Intermediate Crimp



Side Plate
Perforated Sheet

Standard Hybrid Baskets Plain Weave



Outside Dimensions ►

Types ▼	Half Size (1/2 Size)	Middle Size (3/4 Size)	DIN Size	Full Size (1/1 Size)
	Y 396.05B 245 x 255 x 50 mm	Y 296.05B 405 x 255 x 50 mm	Y 196.05B 485 x 255 x 50 mm	Y 197.05B 540 x 255 x 50 mm
	Y 396.06B 245 x 255 x 60 mm	Y 296.06B 405 x 255 x 60 mm	Y 196.06B 485 x 255 x 60 mm	Y 197.06B 540 x 255 x 60 mm
	Y 396.07B 245 x 255 x 70 mm	Y 296.07B 405 x 255 x 70 mm	Y 196.07B 485 x 255 x 70 mm	Y 197.07B 540 x 255 x 70 mm
	Y 396.10B 245 x 255 x 100 mm	Y 296.10B 405 x 255 x 100 mm	Y 196.10B 485 x 255 x 100 mm	Y 197.10B 540 x 255 x 100 mm
	Y 396.05BL 245 x 255 x 50 mm	Y 296.05BL 405 x 255 x 50 mm	Y 196.05BL 485 x 255 x 50 mm	Y 197.05BL 540 x 255 x 50 mm
	Y 396.06BL 245 x 255 x 60 mm	Y 296.06BL 405 x 255 x 60 mm	Y 196.06BL 485 x 255 x 60 mm	Y 197.06BL 540 x 255 x 60 mm
	Y 396.07BL 245 x 255 x 70 mm	Y 296.07BL 405 x 255 x 70 mm	Y 196.07BL 485 x 255 x 70 mm	Y 197.07BL 540 x 255 x 70 mm
	Y 396.10BL 245 x 255 x 100 mm	Y 296.10BL 405 x 255 x 100 mm	Y 196.10BL 485 x 255 x 100 mm	Y 197.10BL 540 x 255 x 100 mm
	Y 392.01 215 x 240 mm	Y 292.01 375 x 240 mm	Y 192.01 455 x 240 mm	Y 193.01 510 x 240 mm

Product Details

Base & Lid
Plain Weave



Side Plate
Perforated Sheet



explore / bahadir.com.tr

Standard Hybrid Baskets Welded Wire Mesh




Length x Width x Height (mm)

Outside Dimensions ►

Types ▼ Lid Base & Lid Lid	Half Size (1/2 Size)	Middle Size (3/4 Size)	DIN Size	Full Size (1/1 Size)
	Y 396.05C 245 x 255 x 50 mm	Y 296.05C 405 x 255 x 50 mm	Y 196.05C 485 x 255 x 50 mm	Y 197.05C 540 x 255 x 50 mm
	Y 396.06C 245 x 255 x 60 mm	Y 296.06C 405 x 255 x 60 mm	Y 196.06C 485 x 255 x 60 mm	Y 197.06C 540 x 255 x 60 mm
	Y 396.07C 245 x 255 x 70 mm	Y 296.07C 405 x 255 x 70 mm	Y 196.07C 485 x 255 x 70 mm	Y 197.07C 540 x 255 x 70 mm
	Y 396.10C 245 x 255 x 100 mm	Y 296.10C 405 x 255 x 100 mm	Y 196.10C 485 x 255 x 100 mm	Y 197.10C 540 x 255 x 100 mm
	Y 396.05CL 245 x 255 x 50 mm	Y 296.05CL 405 x 255 x 50 mm	Y 196.05CL 485 x 255 x 50 mm	Y 197.05CL 540 x 255 x 50 mm
	Y 396.06CL 245 x 255 x 60 mm	Y 296.06CL 405 x 255 x 60 mm	Y 196.06CL 485 x 255 x 60 mm	Y 197.06CL 540 x 255 x 60 mm
	Y 396.07CL 245 x 255 x 70 mm	Y 296.07CL 405 x 255 x 70 mm	Y 196.07CL 485 x 255 x 70 mm	Y 197.07CL 540 x 255 x 70 mm
	Y 396.10CL 245 x 255 x 100 mm	Y 296.10CL 405 x 255 x 100 mm	Y 196.10CL 485 x 255 x 100 mm	Y 197.10CL 540 x 255 x 100 mm
	Y 394.01 215 x 240 mm	Y 294.01 375 x 240 mm	Y 194.01 455 x 240 mm	Y 195.01 510 x 240 mm

Product Details

Base & Lid
Welded Wire Mesh



Side Plate
Perforated Sheet

Specific Hybrid Baskets Plain Weave



Outside Dimensions ▶

Types ▼	Type 1	Type 2 (Dental)	Type 3
	Y 590.01 125 x 110 x 50 mm	Y 398.03 270 x 170 x 30 mm	Y 297.04 340 x 250 x 40 mm
	Y 590.03 125 x 110 x 70 mm	Y 398.05 270 x 170 x 50 mm	Y 298.04 440 x 250 x 40 mm
	Y 590.04 195 x 70 x 40 mm	Y 398.07 270 x 170 x 70 mm	Y 298.08 440 x 250 x 80 mm
	Y 590.05 220 x 125 x 40 mm	Y 398.10 270 x 170 x 100 mm	Y 299.08 460 x 320 x 80 mm
	Y 590.06 265 x 120 x 30 mm		
	Y 590.07 265 x 120 x 40 mm		
	Y 590.08 265 x 120 x 50 mm		
	Y 590.10 265 x 120 x 70 mm		

Product Details

Base & Lid
Plain Weave



Side Plate
Perforated Sheet

Wire Baskets



Contents

	pages
Standard Wire Baskets - Intermediate Crimp	92
Standard Wire Baskets - Plain Weave	93
Standard Wire Baskets - Welded Wire Mesh	94
Optic Wire Baskets	95
Fine Mesh, Mini, Dental Wire Baskets	96
Wide Body & Laundry Wire Baskets	97



BAHADIR

explore / bahadir.com.tr




Intermediate Crimp

Mesh Dimensions: 5x5 mm
Wire Diameter: 1 mm




Plain Weave

Mesh Dimensions: 4x4 mm
Wire Diameter: 1 mm



Welded Wire Mesh

Mesh Dimensions: 6x6mm
Wire Diameter: 1,2 mm



Fine Mesh

Mesh Dimensions: 2x2 mm
Wire Diameter: 1 mm



BAHADIR

explore / bahadir.com.tr

Standard Wire Baskets Intermediate Crimp



Length x Width x Height (mm)

Outside Dimensions ▶



Half Size (1/2 Size)

Y 390.03
245 x 255 x 30 mm

Y 390.05
245 x 255 x 50 mm

Y 390.07
245 x 255 x 70 mm

Y 390.10
245 x 255 x 100 mm

Y 390.01
215 x 240 mm

Middle Size (3/4 Size)

Y 290.03
405 x 255 x 30 mm

Y 290.05
405 x 255 x 50 mm

Y 290.07
405 x 255 x 70 mm

Y 290.10
405 x 255 x 100 mm

Y 290.01
375 x 240 mm

DIN Size

Y 190.03
485 x 255 x 30 mm

Y 190.05
485 x 255 x 50 mm

Y 190.07
485 x 255 x 70 mm

Y 190.10
485 x 255 x 100 mm

Y 190.01
455 x 240 mm

Full Size (1/1 Size)

Y 191.03
540 x 255 x 30 mm

Y 191.05
540 x 255 x 50 mm

Y 191.07
540 x 255 x 70 mm

Y 191.10
540 x 255 x 100 mm

Y 191.01
510 x 240 mm

Product Details



Base & Lid
Intermediate Crimp

Standard Wire Baskets

Plain Weave



Outside Dimensions ▶

Types ▼  Base  Lid	Half Size (1/2 Size)	Middle Size (3/4 Size)	DIN Size	Full Size (1/1 Size)
	Y 392.03 245 x 255 x 30 mm	Y 292.03 405 x 255 x 30 mm	Y 192.03 485 x 255 x 30 mm	Y 193.03 540 x 255 x 30 mm
	Y 392.05 245 x 255 x 50 mm	Y 292.05 405 x 255 x 50 mm	Y 192.05 485 x 255 x 50 mm	Y 193.05 540 x 255 x 50 mm
	Y 392.07 245 x 255 x 70 mm	Y 292.07 405 x 255 x 70 mm	Y 192.07 485 x 255 x 70 mm	Y 193.07 540 x 255 x 70 mm
	Y 392.10 245 x 255 x 100 mm	Y 292.10 405 x 255 x 100 mm	Y 192.10 485 x 255 x 100 mm	Y 193.10 540 x 255 x 100 mm
	Y 392.01 215 x 240 mm	Y 292.01 375 x 240 mm	Y 192.01 455 x 240 mm	Y 193.01 510 x 240 mm

▲ Product Details



Base & Lid
Plain Weave



explore / bahadir.com.tr

Standard Wire Baskets Welded Wire Mesh




Length x Width x Height (mm)

Outside Dimensions ▶


Types ▼  Baklar  Lid	Half Size (1/2 Size)	Middle Size (3/4 Size)	DIN Size	Full Size (1/1 Size)
	Y 394.03 245 x 255 x 30 mm	Y 294.03 405 x 255 x 30 mm	Y 194.03 485 x 255 x 30 mm	Y 195.03 540 x 255 x 30 mm
	Y 394.05 245 x 255 x 50 mm	Y 294.05 405 x 255 x 50 mm	Y 194.05 485 x 255 x 50 mm	Y 195.05 540 x 255 x 50 mm
	Y 394.07 245 x 255 x 70 mm	Y 294.07 405 x 255 x 70 mm	Y 194.07 485 x 255 x 70 mm	Y 195.07 540 x 255 x 70 mm
	Y 394.10 245 x 255 x 100 mm	Y 294.10 405 x 255 x 100 mm	Y 194.10 485 x 255 x 100 mm	Y 195.10 540 x 255 x 100 mm
	Y 394.01 215 x 240 mm	Y 294.01 375 x 240 mm	Y 194.01 455 x 240 mm	Y 195.01 510 x 240 mm

Product Details



Base & Lid
Welded Wire Mesh

Optic Wire Baskets



Outside Dimensions ►



Single Scope

Y 690.02
290 x 80 x 55 mm

Y 690.04
460 x 80 x 55 mm

Y 690.06
670 x 80 x 55 mm

Double Scope

Y 690.03
290 x 80 x 55 mm

Y 690.05
460 x 80 x 55 mm



With Silicone Holders

Outside Dimensions ►



Single Scope

Y 691.02
290 x 80 x 55 mm

Y 691.04
460 x 80 x 55 mm

Double Scope

Y 691.03
290 x 80 x 55 mm

Y 691.05
460 x 80 x 55 mm



With Silicone Holders


Wire Range - 13,5 mm

Optic wire baskets are made from stainless steel.
Diameters: 2 mm - 2,5 mm - 3 mm - 4 mm.



explore / bahadir.com.tr

Standard Wire Baskets Mini Model & Dental Model



Length x Width x Height (mm)

Outside Dimensions ▶



Mini Model (Intermediate Crimp)

Y 490.05
265 x 120 x 45 mm



Mini Model (Plain Weave)

Y 492.05
265 x 120 x 45 mm



Outside Dimensions ▶




Dental Model (Intermediate Crimp)

Y 590.12
270 x 170 x 45 mm




Dental Model (Plain Weave)

Y 590.13
270 x 170 x 45 mm



/ Fine Mesh Baskets

Outside Dimensions



Y 591.01
80 x 80 x 40 mm

Y 591.03
145 x 150 x 50 mm

Y 591.07
220 x 140 x 50 mm

Y 591.08
275 x 175 x 35 mm


Y 591.09
290 x 170 x 60 mm

Product Details



© Bahadır 2020

Wide Body Wire Baskets



Outside Dimensions

Types

▼



Wide Body
(Intermediate Crimp)

Y 190.55
470 x 337 x 51 mm

Y 190.60
533 x 337 x 101 mm



Outside Dimensions



/ Surgical Laundry Baskets

Types

▼



Y 693.02
605 x 405 x 150 mm

Y 693.04
605 x 405 x 300 mm

pages

Container Accessories

100-101

Wire Basket & Perforated Tray Accessories

102-106

CONTENTS




ACCESSORIES




Filter Retainer

Y 808.01




Filter Retainer
(For M Series)

For Paper Filter - Y 808.00




Filter Dispenser
(For Round Filters)

Y 801.00



Filter Retainer
(Aluminum)

Y 805.00



Filter Retainer
(Aluminum)

Y 805.10



Filter Retainer
(Stainless Steel)

Y 805.11



Filter Retainer
(Stainless Steel)

Y 805.20



Bio Barrier
(Cassette Filter)

Y 806.00B



Valve
(Cassette Filter)

Y 806.00V



PTFE
(Cassette Filter)

Y 806.00P

Identification Labels

Gynaecology	Engraved
Spinal	Y 820.01 ●
ENT	Y 820.02 ●
General	Y 820.03 ●
Neurosurgery	Y 820.04 ●
Micro	Y 820.05 ●
	Y 820.06 ●

Not Engraved
Y 820.11 ●
Y 820.12 ●
Y 820.13 ●
Y 820.14 ●
Y 820.15 ●
Y 820.16 ●

Gynaecology	Engraved
Spinal	Y 820.01M ●
ENT	Y 820.02M ●
General	Y 820.03M ●
Neurosurgery	Y 820.04M ●
Micro	Y 820.05M ●
	Y 820.06M ●

Not Engraved
Y 820.11M ●
Y 820.12M ●
Y 820.13M ●
Y 820.14M ●
Y 820.15M ●
Y 820.16M ●



Paper Filter
(For 1/1, 3/4, 1/2 Size)

100 pcs - Y 800.00



Paper Filter
(For Mini Container)

100 pcs - Y 802.00



Paper Filter
(For Dental Container)

100 pcs - Y 804.00



Paper Filter
(For Micro Dental Container)

100 pcs - Y 803.00



PTFE Filter
(For 1/1, 3/4, 1/2 Size)

Y 800.10



PTFE Filter
(For Mini Container)

Y 802.10



PTFE Filter
(For Dental Container)

Y 802.11




PTFE Filter
(For Micro Dental Container)

100 pcs - Y 803.10



Textile Filter
(For 1/1, 3/4, 1/2 Size)

Y 800.01



Textile Filter
(For Mini & Dental Container)

Y 802.01




Textile Filter
(For Micro Dental Container)

Y 803.01



Label Dispenser
(For Identification Labels)

Y 829.00



Indicators
(Paper)

100 pcs - Y 810.00

100 pcs - Y 810.01

100 pcs - Y 812.00



Security Seal

100 pcs - Y 819.00



Security Seal
(Service)

100 pcs - Y 819.01



Silicone Mat
(Rectangle - Unperforated)

Green - Y 850.00

Blue - Y 850.01

270 x 125 mm



Silicone Mat
(Rectangle - Perforated)

Green - Y 850.10

Blue - Y 850.11

270 x 125 mm



Silicone Mat
(Square - Unperforated)

Green - Y 851.00

Blue - Y 851.01

245 x 245 mm



Silicone Mat
(Square - Perforated)

Green - Y 851.10

Blue - Y 851.11

245 x 245 mm



Silicone Mat
(Green - Perforated)

225 x 235 mm - Y 852.00

370 x 235 mm - Y 852.10

450 x 235 mm - Y 852.20

500 x 235 mm - Y 852.30




Silicone Mat
(Blue - Perforated)

225 x 235 mm - Y 852.01

370 x 235 mm - Y 852.11

450 x 235 mm - Y 852.21

500 x 235 mm - Y 852.31





It can be cut, if desired.

Silicone Mat
(Elevated Nubs - Perforated)

Y 853.00

545 x 255 mm





Types

For Optic
Wire BasketFor Optic
Perforated TraySilicone Holder
Ø 12 mmY 856.01
37 x 29 x 13 mmSilicone Holder
Ø 14 mmY 856.02
37 x 29 x 13 mmY 856.05
30 x 29 x 13 mmY 856.06
30 x 29 x 13 mm

Silicone Holder for Endoscopic Instruments

Silicone holders provide a safe sterilization process
for endoscopic instruments.

Endoscopic Rack

485 x 235 x 160 mm
(Y 855.01 x 6 pcs.)
(Y 855.02 x 24 pcs.)Silicone Holder
(For Endoscopy Racks)

Ø 10 mm - Y 855.01

Ø 5 mm - Y 855.02

Complete - Y 843.06

Silicone Holder
(For Optic Wire Baskets)

40 x 17 x 6,5 mm - Y 858.03

70 x 17 x 6,5 mm - Y 858.04

Complete - Y 858.00

Silicone Holder
(For Optic Wire Baskets)

54 x 17 x 6,5 mm - Y 859.04

68 x 17 x 6,5 mm - Y 859.06

80,5 x 17 x 6,5 mm - Y 859.07

Complete - Y 859.00

Silicone Holder
(For Adjustable Racks)

226 x 17 mm - Y 857.01

226 x 17 mm - Y 857.03

Complete - Y 843.04

230 x 220 x 40 mm
(Y 857.01 x 1 pcs.)
(Y 857.03 x 2 pcs.)

Adjustable Rack





Silicone Rack
(Type 1)

Blue - Y 854.00
Yellow - Y 854.00S



Blue (Complete) - Y 854.70

(8 x 26 x 40 mm)
(Y 854.00 x 4 pcs.)
(Y 854.60P x 1 pcs.)



Yellow (Complete) - Y 854.70S

(8 x 26 x 40 mm)
(Y 854.00S x 4 pcs.)
(Y 854.60P x 1 pcs.)



Silicone Rack
(Type 2)

Blue - Y 854.03
Yellow - Y 854.03S



Blue (Complete) - Y 854.73

(8 x 26 x 40 mm)
(Y 854.03 x 4 pcs.)
(Y 854.60P x 1 pcs.)



Yellow (Complete) - Y 854.73S

(8 x 26 x 40 mm)
(Y 854.03S x 4 pcs.)
(Y 854.60P x 1 pcs.)



Silicone Rack
(Type 3)

Blue - Y 854.05
Yellow - Y 854.05S



Blue (Complete) - Y 854.75

(8 x 26 x 40 mm)
(Y 854.05 x 4 pcs.)
(Y 854.60P x 1 pcs.)



Yellow (Complete) - Y 854.75S

(8 x 26 x 40 mm)
(Y 854.05S x 4 pcs.)
(Y 854.60P x 1 pcs.)



Polymer Rail
(16 x 100 x 9 mm)

Y 854.63P




Polymer Clip
(16 x 100 x 9 mm)

Y 854.61P




Complete - Y 854.60P

(Y 854.63P x 1 pcs.)
(Y 854.61P x 4 pcs.)



Polymer rails and polymer clips can be used for all types of wire baskets and perforated trays.



Holding Pin
(For Wire Baskets)

Y 830.00



Holding Pin
(For Partition Sheets)

Y 830.01



Partition Sheets
(Stainless Steel)

50x20 mm - Y 831.05

130x20 mm - Y 831.13

225x20 mm - Y 831.22

460x20 mm - Y 831.46



Silicone Divider

Y 840.00



Sliding Rails
(With Screws)

230 mm - Y 841.00


470 mm - Y 841.01



Limiting Bars
(With Screws)

230 mm - Y 842.00

470 mm - Y 842.01



Spiral Holder
(Single)

Y 845.00



Spiral Holder
(Double)


Y 847.00



Fixation Clamps

50x20x12 mm - Y 849.00

65x20x20 mm - Y 849.01



Holding Clamps

24x28x25 mm - Y 849.20

Ø 7-12 mm

24x28x29 mm - Y 849.21

Ø 12-16 mm

24x28x31 mm - Y 849.22

Ø 16-26 mm

24x28x40 mm - Y 849.23

Ø 26-36 mm

24x28x44 mm - Y 849.24

Ø 36-45 mm



Small tag for wire baskets



Tag with slot for RFID Placement



Small Tag
(For Wire Baskets)

40 x 20 mm - Y 837.00




Large Tag
(For Wire Baskets)

50 x 30 mm - Y 837.01



Slotted Tag
(For Wire Baskets)

52 x 18 mm - Y 838.00



Silicone Dividers
(With Rail)

30 mm - Y 854.76

115 mm - Y 854.78

235 mm - Y 854.80



Silicone Divider

235 mm - Y 854.30



Stainless Steel Dividers

110 x 38 mm - Y 842.05

235 x 38 mm - Y 842.06



Surgical Drapes
(Cotton)

60 x 40 mm - Y 095.01

65 x 65 mm - Y 095.02

85 x 85 mm - Y 095.03

90 x 60 mm - Y 095.04

90 x 80 mm - Y 095.05

130 x 45 mm - Y 095.06

130 x 60 mm - Y 095.07

140 x 100 mm - Y 095.08

185 x 45 mm - Y 095.09

185 x 60 mm - Y 095.10

145 x 140 mm - Y 095.17



Drape Frames

545 x 255 mm - Y 195.00

405 x 255 mm - Y 295.00

255 x 255 mm - Y 395.00

280 x 125 mm - Y 495.00

310 x 190 mm - Y 595.00



Surgical Trolley

1157x706x1405 mm - Y 890.00