

## Skeletal traction, simplified.



### IMPROVE PATIENT'S CARE

Alleviate the patient's pain faster while maintaining sterility.



### NO MORE REPROCESSING

All the tools you need are pre-sterilized to save you time and money.




### COMPLIMENTARY SURGICAL POWER KIT

No need to borrow OR power, we provide a surgical drill kit, at no cost.




### CONSISTENCY ON THE PROCEDURE

We provide the surgical drill and training to standardize the procedure.



# Skeletal Traction Procedure Kit

Everything you need for skeletal traction, in one package



**Product Name:** Skeletal Traction Procedure Kit **Catalogue #:** TP-003-06

**Product Description:** 1 SteriTrak® pin, 2mm, 12.2 inch k-wire, 1 set of pin tip covers, 1 SteriGo Single-use Cover, 4 Blue Towels, 4 Sticky drapes, 1 Roll of Kerlix, 1 20cc Syringe, 2 18ga x 1.5 Needle, 2 22ga x 1.5 Needle, 1 11 Blade, 1 Skin marking pen, 2 4"x4" gauze pads, 2 Xeroform pads, 1 Traction rope, 2 10.5ml of Chloraprep w/ Applicator



**Product Name:** Complimentary Surgical Power Kit **Catalogue #:** ST-001

**Product Description:** 1 modified FDA cleared DeWalt drill, 2 batteries, 1 charger, 1 pin snip, and 1 drill carry bag.

## Proud partner of Stanley Black & Decker since 2017

Arbutus Medical developed a custom surgical drill, powered by DEWALT®, to offer you the market's most reliable pin driver for skeletal traction. We understand that power is not always available outside the OR, that's why, we provide a compact power kit specifically designed for skeletal traction at no additional cost.

