


764-550 Rod and Pin Cutter 22" Length Maximum Capacity-6.3mm



Indications for use: The Amerisurg Rod Cutter is designed to be used as an invasive surgical instrument for cutting spinal rods and pins, in situ.

Standard Cleaning and Sterilization instructions apply to ensure instrument is properly cleaned and sterilized based on hospital policy and minimum industry standards.

The Amerisurg Rod Cutter arms are detachable for easier cleaning and sterilization.

The tip of the jaws may be customized based on surgeon need, i.e. more tapered tip, more rounded tip etc. Consult your salesman for further details