



American Surgical

Specialties Company



A CATALOG
OFFERING
SPECIALTY PRODUCTS
FOR THE
OPERATING ROOM
&
STERILE PROCESSING



 **AmeriSurg**

WOMEN OWNED & OPERATED

www.americansurg.com



Designed For You By You

Customer Satisfaction has been and will continue to be our top priority. To those nurses who helped with the development of this catalog, *we thank you!* Your comments are very much appreciated: “Easy One-Stop-Shop”, “Great things come in small catalogs”, “Just what is needed for me in one book.”

TABLE OF CONTENTS

Tip Protectors	Page 1 - 4
Clamp Covers	Page 5
Suture Booties	Page 5
ID Tape	Page 6 - 7
Vessel Loops	Page 8
Instrument Stringers	Page 8
Light Handle Covers	Page 9
Trays/Pans	Page 10 - 13
Cleaning Brushes/Brush Kits	Page 14 - 16
MonoFoil	Page 17
Instrument Care Soap/Lubrication	Page 18 - 19
Alton Dean Pump	Page 20
Trump-It Irrigation	Page 20
Insufflation Tubing & Accessories	Page 21 - 22
Scope Warmer	Page 22
Endoscopic Fiberoptic Light Cables	Page 23
Laparoscopic Cutting Loops	Page 24
Urological Holmium Fibers	Page 24
Laser Guard	Page 25 - 27
OB/Gyn Uterine Manipulator & Dilators	Page 28 - 29
Lletz/LEEP Electrodes & Accessories	Page 30 - 32
Bovie Electrodes & Cauteries	Page 33 - 34
EYE: Ophthalmic Burrs	Page 35
Nerve Locators	Page 35
Lasiks UltraCell Products	Page 36
Cataract/Refractive Products	Page 37
UltraCell ENT	Page 38 - 40
EYE Knives & Blades	Page 41 - 43
Trephines	Page 43
Cataract Cannulas	Page 44 - 47
ENT Knives & Blades	Page 48 - 50
Surgical Marking Pens	Page 51
Drapes: C-Arm	Page 52 - 53
Drapes: Urology/Ultrasound	Page 54
Drapes: Microscope Drapes	Page 55
Shielding-Lead Aprons	Page 56
Shielding-Lead Accessories	Page 57

Standard Instrument Tip Protectors

Standard Instrument Tip Protectors are available solid, tinted or clear and vented or non-vented. They are gas or steam permeable and can be used in flash or conventional autoclave processes. All are available in 100 per package. Latex-free and single-use.



1.6mm x 19mm

Solid Vented	091011
Solid Non-Vented	091001
Tinted Vented	093011
Tinted Non-Vented	093001
Clear Vented	092011
Clear Non-Vented	092001



2mm x 19mm

Solid Vented	091012
Solid Non-Vented	091002
Tinted Vented	093012
Tinted Non-Vented	093002
Clear Vented	092012
Clear Non-Vented	092002



2.8mm x 19mm

Solid Vented	091013
Solid Non-Vented	091003
Tinted Vented	093013
Tinted Non-Vented	093003
Clear Vented	092013
Clear Non-Vented	092003



3.2mm x 25.4mm

Solid Vented	091014
Solid Non-Vented	091004
Tinted Vented	093014
Tinted Non-Vented	093004
Clear Vented	092014
Clear Non-Vented	092004



5mm x 25mm

Solid Vented	091015
Solid Non-Vented	091005
Tinted Vented	093015
Tinted Non-Vented	093005
Clear Vented	092015
Clear Non-Vented	092005



10mm x 19mm

Solid Vented	091016
Solid Non-Vented	091006
Tinted Vented	093016
Tinted Non-Vented	093006
Clear Vented	092016
Clear Non-Vented	092006



2mm x 9mm x 25mm

Solid Vented	091017
Solid Non-Vented	091007
Tinted Vented	093017
Tinted Non-Vented	093007
Clear Vented	092017
Clear Non-Vented	092007



2mm x 16mm x 25mm

Solid Vented	091018
Solid Non-Vented	091008
Tinted Vented	093018
Tinted Non-Vented	093008
Clear Vented	092018
Clear Non-Vented	092008



3mm x 25mm x 25mm

Solid Vented	091019
Solid Non-Vented	091009
Tinted Vented	093019
Tinted Non-Vented	093009
Clear Vented	092019
Clear Non-Vented	092009

Assorted — Solid

Solid Vented	091020
Solid Non-Vented	091010

Assorted — Tinted

Tinted Vented	093020
Tinted Non-Vented	093010
Tinted Vented Micro	093000

Assorted — Clear

Clear Vented	092020
Clear Non-Vented	092010

DuoGuards™

Ensure that your instruments stay in the proper open position during processing with DuoGuards™. Available vented and non-vented, DuoGuards™ are tinted for easy viewing of instrument tips. Latex-free and single-use.



Tinted
50 Per Bag

Part No. Vented	Part No. Non-Vented	Color	mm
093031	093021	White	6.4 x 15
093032	093022	Blue	9.5 x 20
093033	093023	Green	12.7 x 20
093034	093024	Red	19 x 25.5
093035	093025	Yellow	15 x 25
093037	093027	Orange	25 x 25
093040	093030	Assorted	Assorted

Low Temperature Web Guards

Another exclusive, these unique guards with a flow through, open weave are designed to be used exclusively with low heat, plasma sterilization processes, such as the Sterrad™ Sterilization System. They are great for protecting endoscopes and their accessories. Latex-free and single-use.



*Sterrad is a Trademark of Johnson & Johnson Medical, Inc.

Part No.	Color	mm	Quantity
099001	Orange Mesh	10 x 76	24
099002	Green Mesh	15 x 76	24
099003	Blue Mesh	20 x 76	24
099004	Red Mesh	40 x 76	24
099005	Black Mesh	60 x 76	16
099006	White Mesh	70 x 76	16
099007	Assorted Mesh	Assorted	20

InnerVent™ Instrument Tip Protectors

InnerVent™ Instrument Tip Protectors offer the performance of vented tip protectors with the safety of non-vented tip protectors. The raised inner nubs allow improved flow through of steam and gas while preventing the tip protectors from "vacuuming" onto the instrument during drying. Latex-free and single-use.



Tinted



Solid/Radiopaque



Clear

Solid, Clear, Tinted
100 Per Bag

Part No.	Color	mm
091117	Orange, Solid	2 x 9 x 25
091118	Brown, Solid	2 x 16 x 25
091119	Black, Solid	3 x 25 x 25
092117	Clear	2 x 9 x 25
092118	Clear	2 x 16 x 25
092119	Clear	3 x 25 x 25
093117	Orange, Tinted	2 x 9 x 25
093118	Brown, Tinted	2 x 16 x 25
093119	Black, Tinted	3 x 25 x 25

Gelpi Retractor Tip Protectors

Protect expensive Gelpi retractors with exclusive Gelpi Retractor Tip Protectors. Available in 3 sizes to fit the most commonly used Gelpi retractors, these tip protectors are gas and steam sterilization compatible. Latex-free and single-use.



Part No.	Color	mm	Quantity Per Bag
096109	Clear	38 Length	55
096129	Clear	51 Length	45
096139	Clear	76 Length	40

Weitlaner Instrument Tip Protectors

The Weitlaner Retractor Tip Protector is specifically designed for use with up to a 4" Weitlaner retractor and offers easy one-handed application. Protect your expensive Weitlaner instruments and reduce reprocessing costs associated with compromised wrap. Latex-free and single-use.



Vented, Solid or Clear
50 Per Bag

Part No.	Size	Color
096409	Small	Magenta
096410	Small	Clear
096419	Medium	Green
096420	Medium	Clear

Laparoscopic MIS Tip Protectors

Venting of these specialty tip protectors provides extra flow through to Trocar probes and Verres needles during high level disinfection or sterilization. Latex-free and single-use.



Vented and Clear

Part No.	Applications	Quantity Per Bag
097011	19mm, Use for Verres Needle	50
097012	20mm 5mm, Use for Trocar Probe	25
097013	43mm 5mm, Use for Trocar Probe	20
097014	20mm 10mm, Use for Trocar Probe	25
097015	57mm 10mm, Use for Trocar Probe	20

Endoscope Guards

Safeguard the eyepieces of delicate endoscopic equipment during storage with these soft, conformable guards. Latex-free and single-use.



Blue, Non-Vented
6 Per Bag

Part No.	Description
095016	Small Eye Piece

Round Tip Protectors

Designed to protect valuable instruments ranging in sizes from the smallest Verres Needle to the largest endoscope during handling and storage. Tinted for easy visual inspection and vented for improved steam, gas, or disinfectant flow. Latex-free and single-use.



Vented and Tinted

Part No.	Part No.	Color	Fits Tip Sizes mm	Quantity Per Bag
096001	096021	White	2-3	70
096002	096022	Blue	3-4	70
096003	096023	Green	4-5	70
096004	096024	Red	5-6	70
096005	096025	Yellow	6-7	70
096006	096026	Maroon	7-8	70
096007	096027	Orange	8-9	70
096008	096028	Brown	9-10	70
096009	096029	Black	10-11	70
096013	096033	Green	15-16	50
096015	096035	Yellow	19-20	50
096017		Orange	21-22	50

Instrument Tip Protectors

Instrument Tip Protectors are designed to protect valuable instruments from breaking and chipping during sterilization and handling, help eliminate reprocessing costs caused by torn pouches or wrap, and protect staff from sharp sticks during handling.

Osteotome Tip Protectors

Save Osteotome instrument edges from chipping with these tip protectors. Available tinted or clear and vented or non-vented to comply with user preferences. All types allow visual instrument inspection and are gas and steam permeable. Latex-free and single-use. Available in clear or blue.



Clear
50 Per Bag

Part No. Vented	Part No. Non-Vented	mm
094031	094021	6.4 x 19
094032	094022	9.5 x 25
094033	094023	12.7 x 25
094034	094024	19 x 25.4
094035	094025	25 x 25
094036	094026	32 x 32
094037	094027	38 x 38
094038	094028	44.5 x 32
094040	094030	Assorted




Blue
50 Per Bag

Part No. Vented	Part No. Non-Vented	mm
094011	094001	6.4 x 19
094012	094002	9.5 x 25
094013	094003	12.7 x 25
094014	094004	19 x 25.4
094015	094005	25 x 25
094016	094006	32 x 32
094017	094007	38 x 38
094018	094008	44.5 x 32
094020	094010	Assorted

Standard Tip Protector Sizing Guide

Use this guide to determine the standard tip protector that will best fit your instrument.



Clamp Covers


Clamp Covers assist in the atraumatic clamping of delicate vessels during surgery. They are available in radiopaque silicone, non-radiopaque fabric, and radiopaque fabric. We offer a variety of sizes to provide outstanding convenience and flexibility of use. These are packaged sterile or non-sterile, intended for single-use, and latex-free.



Part No.	Color	Material	Radiopaque	Length mm	Length Inches	Quantity Per Box
Sterile						
061011	Red	Silicone	Yes	2.5 x 130	0.10 x 5	10 Pkgs of 2
061012	White	Silicone	Yes	5.1 x 130	0.20 x 5	10 Pkgs of 2
061013	Blue	Silicone	Yes	7.6 x 130	0.30 x 5	10 Pkgs of 2
071001	White	Fabric	No	3.2 x 130	0.125 x 5	10 Pkgs of 4
071002	White	Fabric	No	9.5 x 130	0.375 x 5	10 Pkgs of 4
072001	White	Fabric	Yes	3.2 x 150	0.125 x 6	10 Pkgs of 4
072002	White	Fabric	Yes	9.5 x 150	0.375 x 6	10 Pkgs of 4
072004	White	Fabric	Yes	6.3 x 150	0.25 x 6	10 Pkgs of 4
Non-Sterile						
062011	Red	Silicone	Yes	2.5 x 130	0.10 x 5	20
062012	White	Silicone	Yes	5.1 x 130	0.20 x 5	20
062013	Blue	Silicone	Yes	7.6 x 130	0.30 x 5	20

Suture Aid Booties

Suture Aid Booties protect, tag, and help locate sutures. Our vinyl suture booties are available in two-, three-, or five-pair configurations and are positioned in an adhesive-backed foam block for quick and easy usage. They are radiopaque, latex-free, and are packaged sterile for single-use. Choose from a variety of sizes and colors as listed below.



Part No.	Description	Quantity Per Tray	Quantity Per Box
051003	Standard Yellow	5 Pair	5 Trays
051011	Standard Red	5 Pair	5 Trays
051012	Standard Blue	5 Pair	5 Trays
051014	Standard White	5 Pair	5 Trays
051015	Standard Asst.	5 Pair	5 Trays
051004	Large Pink	5 Pair	5 Trays
053003	Standard Yellow	3 Pair	5 Trays
052004	Large Pink	2 Pair	5 Trays
052003	Standard Yellow	2 Pair	5 Trays

Instrument Identification Tape

Tape makes organizing and tracking valuable instruments easy! Tape is flexible, adheres securely to instruments, and will not chip, flake, peel, or crack — and with 75+ colors, symbols, and character options, the color-coding options are nearly endless. If a change is needed, simply remove the ColorTrack I.D. Tape and replace with your new color choice. Latex-free.

Sheet Tape

Each 8½" x 11" sheet has 374 total inches of tape and contains:

- 44 each of 2 inch strips
- 44 each of 1½ inch strips
- 220 each of 1 inch strips

Solid Color

¼" Height / Various widths

Part No.	Description
152001	White
152002	Royal Blue
152003	Kelly Green
152004	Red
152005	Yellow
152006	Purple
152007	Orange
152008	Brown
152009	Black
152010	Gray
152011	Lavender
152012	Beige
152013	Light Blue
152014	Lime Green
152015	Pink
152017	Plum
152018	Mint Green

Dots

¼" Height / Various widths

Part No.	Description
158040	Dots, Fuchsia
158041	Dots, Medical Blue
158042	Dots, Emerald Green
158043	Dots, Violet

Stars

¼" Height / Various widths

Part No.	Description
158030	Stars, Fuchsia
158031	Stars, Medical Blue
158032	Stars, Emerald Green
158033	Stars, Violet

White Stripe

¼" Height / Various widths

Part No.	Description
154002	White/Royal Blue
154003	White/Kelly Green
154004	White/Red
154005	White/Yellow
154007	White/Orange
154008	White/Brown
154011	White/Lavender
154013	White/Light Blue
154014	White/Lime Green
154015	White/Pink
154017	White/Plum

Diagonal Stripes

¼" Height / Various widths

Part No.	Description
159001	Diagonal, Purple/Black
159002	Diagonal, Orange/Black
159003	Diagonal, Blue/Black
159004	Diagonal, Rubine/Black
159005	Diagonal, Yellow/Black
159006	Diagonal, Green/Black
159007	Diagonal, White/Black

Tri-Color Stripe

¼" Height / Various widths

Part No.	Description
155051	Red/White/Blue
155052	Red/White/Black
155053	Red/Yellow/Black
155054	Red/Black/Green

Black Stripe

¼" Height / Various widths

Part No.	Description
153001	Black/White
153003	Black/Kelly Green
153004	Black/Red
153005	Black/Yellow
153007	Black/Orange
153010	Black/Gray
153011	Black/Lavender
153012	Black/Beige
153013	Black/Light Blue
153014	Black/Lime Green
153015	Black/Pink
153017	Black/Plum

Triangles

¼" Height / Various widths

Part No.	Description
158050	Triangles, Fuchsia
158051	Triangles, Medical Blue
158052	Triangles, Emerald Green
158053	Triangles, Violet

Numbers

¼" Height / Various widths

Each sheet contains numbers 0–8

Part No.	Description
158020	Numbers, Fuchsia
158021	Numbers, Medical Blue
158022	Numbers, Emerald Green
158023	Numbers, Violet

Sheet Tape continued

Assorted Symbols

¼" Height / Various widths

Part No.	Description
158060	Red Hearts on White
158070	Smiley Face, Yellow/Black
158080	Alphabet, Fuschia/Wt, A-E
158061	White Hearts on Red

4-Color

¼" Height / Various widths

Part No.	Description
152031	Red/Yellow/Purple/Orange
152032	Blue/White/Black/Green

Specialty Stripe

¼" Height / Various widths

Part No.	Description
156003	Blue/Green
156004	Blue/Red
157002	Yellow/Blue
157003	Yellow/Green
157004	Yellow/Red

Code to Specialty

¼" Height / Various widths

Part No.	Description
ENT 158090	Mint/White ENT
ER ER 158091	Purple/White ER
GEN 158092	Royal Blue/White GEN
OB/GYN 158093	Lt Blue/White OB/GYN
ORTHO 158094	Kelly Green/White ORTHO
NEURO 158095	Orange/White NEURO
CARDIAC 158096	Red/ White CARDIAC

Roll Tape

Description	1/8" x 250" Part No.	1/4" x 250" Part No.	3/8" x 250" Part No.
White	AMT-01-01	AMT-02-01	AMT-03-01
Royal Blue	AMT-01-02	AMT-02-02	AMT-03-02
Kelly Green	AMT-01-03	AMT-02-03	AMT-03-03
Red	AMT-01-04	AMT-02-04	AMT-03-04
Yellow	AMT-01-05	AMT-02-05	AMT-03-05
Purple	AMT-01-06	AMT-02-06	AMT-03-06
Orange	AMT-01-07	AMT-02-07	AMT-03-07
Brown	AMT-01-08	AMT-02-08	AMT-03-08
Black	AMT-01-09	AMT-02-09	AMT-03-09



Color Coded Identification Locking Tags

Tamper-resistant. Easily removed without scissors. Non-sterile, Autoclavable, Latex-free, Versatile, Durable and Economical.

1" Length. 100 per package



Description	Part No.
White	AMLT-100
Blue	AMLT-200
Green	AMLT-300
Red	AMLT-400
Yellow	AMLT-500
Black	AMLT-900
Pink	AMLT-110
Gold	AMLT-112
Purple	AMLT-116
Grey	AMLT-600
Orange	AMLT-700
Brown	AMLT-800

Vessel Loops

Vessel loops are used in occluding, retracting, and identifying arteries, veins, tendons and nerves in surgical procedures. Our vessel loops are extruded from medical-grade silicone. All vessel loops are latex-free, radiopaque, and intended for single-use



2 Loops / Pkg; 10 Pkgs / Box Double Pouched

Part No.	Color	mm Width x Length	Inches Width x Length
011001	Red Mini	0.83 x 406	0.03 x 16
011002	Blue Mini	0.83 x 406	0.03 x 16
011003	Yellow Mini	0.83 x 406	0.03 x 16
011004	White Mini	0.83 x 406	0.03 x 16
011011	Red Maxi	1.0 x 406	0.04 x 16
011012	Blue Maxi	1.0 x 406	0.04 x 16
011013	Yellow Maxi	1.0 x 406	0.04 x 16
011014	White Maxi	1.0 x 406	0.04 x 16
011021	White Supermaxi	1.3 x 559	0.05 x 22
011022	Blue Supermaxi	1.3 x 559	0.05 x 22
011032	Blue Micro	0.41 x 559	0.019 x 22
011301	Red Mini	0.83 x 762	0.03 x 30
011302	Blue Mini	0.83 x 762	0.03 x 30
011303	Yellow Mini	0.83 x 762	0.03 x 30
011304	White Mini	0.83 x 762	0.03 x 30
011311	Red Maxi	1.0 x 762	0.04 x 30
011312	Blue Maxi	1.0 x 762	0.04 x 30
011313	Yellow Maxi	1.0 x 762	0.04 x 30
011314	White Maxi	1.0 x 762	0.04 x 30

2 Loops / Pkg; 10 Pkgs / Box Carded and Pouched

Part No.	Color	mm Width x Length	Inches Width x Length
016001	Red Mini	0.83 x 457	0.03 x 18
016002	Blue Mini	0.83 x 457	0.03 x 18
016003	Yellow Mini	0.83 x 457	0.03 x 18
016004	White Mini	0.83 x 457	0.03 x 18
016011	Red Maxi	1.0 x 457	0.04 x 18
016012	Blue Maxi	1.0 x 457	0.04 x 18
016013	Yellow Maxi	1.0 x 457	0.04 x 18
016014	White Maxi	1.0 x 457	0.04 x 18

Light Handle Products

Surgical Light Handle Covers

- Increases productivity, reduces setup time
- Provides effective sterile barrier
- Light handle adapters fit all brands
- Available sterile or non-sterile

Flexible Handle

Sterile



Metal handle not included but available for Flexible & Semi-Rigid Handles

Part No.	Description
30-0511	1/Pkg. 160/Case
30-0512	2/Pkg. 80/Case
30-0513	3/Pkg. 64/Case

Semi-Rigid Handle

Sterile



Part No.	Description
30-0711	1/Pkg. 40/Case
30-0712	2/Pkg. 20/Case

Rigid Handle

Sterile



Adapter not included but available

Part No.	Description
30-0881	1/Pkg. 64/Case
30-0882	2/Pkg. 40/Case

Instrument Stringers

- Promotes safer, quicker instrument handling, better processing and decontamination
- Sets sterile field quicker
- Many sizes and styles available



Spring Loaded

Part No.	Inches L x W
10-1000	6"L x 3"W
10-1001	8"L x 3"W
10-1002	10"L x 3"W
10-1003	12"L x 3"W
10-1004	14"L x 3"W
10-1005	16"L x 3"W



Pull to Open Decontamination

Part No.	Inches L x W
10-1021	6"L x 5.5"W
10-1022	8"L x 5.5"W
10-1023	10"L x 5.5"W
10-1024	12"L x 5.5"W



Center Locking

Part No.	Inches L x W
10-1006	6"L x 2.5"W
10-1007	8"L x 2.5"W
10-1008	10"L x 2.5"W
10-1009	12"L x 2.5"W
10-1010	14"L x 2.5"W
10-1011	16"L x 2.5"W



Open Ended U-Shaped

Part No.	Inches L x W
10-1025	6"L x 2.75"W
10-1026	8"L x 2.75"W
10-1027	10"L x 2.75"W
10-1028	12"L x 2.75"W
10-1029	14"L x 2.75"W
10-1030	16"L x 2.75"W



Pull to Open

Part No.	Inches L x W
10-1012	4"L x 2.5"W
10-1013	6"L x 2.5"W
10-1014	8"L x 2.5"W
10-1015	10"L x 2.5"W
10-1016	12"L x 2.5"W
10-1017	14"L x 2.5"W
10-1018	16"L x 2.5"W
10-1019	18"L x 2.5"W
10-1020	20"L x 2.5"W



Open Ended U-Shaped

Part No.	Inches L x W
10-1031	6"L x 5"W
10-1032	8"L x 5"W
10-1033	10"L x 5"W
10-1034	12"L x 5"W
10-1035	14"L x 5"W
10-1036	16"L x 5"W

INSTRUMENT TRAYS & FLAT/SOLID COVERS

- Wide assortment of heavy duty, seamless trays have covered corners for easy sterilization
- Covers have recessed strap handles for ease in handling and storing

	DIMENSIONS	CASE	CASE
	L x W x H	LOT	WEIGHT
AM43	5 $\frac{1}{2}$ " x 4" x 2 $\frac{1}{2}$ "	12	5.0 lbs
AM43-2	Cover—Fits 43 Tray	12	3.5 lbs
AM63	6 $\frac{1}{2}$ " x 6 $\frac{1}{2}$ " x 2 $\frac{1}{4}$ "	12	7.3 lbs
AM63-2	Cover—Fits 63 Tray	12	5.1 lbs
AM952	8 $\frac{1}{2}$ " x 5" x 2"	12	7.3 lbs
AM953	8 $\frac{1}{2}$ " x 5" x 3"	12	8.6 lbs
AM952-2	Cover—Fits 952 & 953 Trays	12	5.8 lbs
AM1002	10" x 6 $\frac{1}{2}$ " x 2"	12	10.3 lbs
AM1002-2	Cover—Fits 1002 Tray	12	8.5 lbs
AM1202	12 $\frac{1}{2}$ " x 7 $\frac{1}{2}$ " x 2 $\frac{1}{4}$ "	12	16.0 lbs
AM1202-2	Cover—Fits 1202 Tray	12	11.3 lbs
AME1082	10" x 8 $\frac{1}{2}$ " x 2 $\frac{1}{2}$ "	6	7.2 lbs
AME1082-2	Cover—Fits 1082 Tray	6	6.7 lbs
AME1650	16 $\frac{1}{2}$ " x 10" x 2 $\frac{1}{2}$ "	6	12.5 lbs
AME1654	16 $\frac{1}{2}$ " x 10" x 4"	4	9.7 lbs
AM1650-2	Cover—Fits Above Two Trays	6	9.3 lbs
AME902	6 $\frac{1}{2}$ " x 4 $\frac{1}{2}$ " x 2"	12	6.0 lbs
AME904	6 $\frac{1}{2}$ " x 4 $\frac{1}{2}$ " x 4"	12	8.1 lbs
AM902-2	Cover—Fits Above Two Trays	12	2.9 lbs
AME2302	13 $\frac{1}{2}$ " x 12 $\frac{1}{2}$ " x 2 $\frac{1}{2}$ "	6	14.3 lbs
AME2304	13 $\frac{1}{2}$ " x 12 $\frac{1}{2}$ " x 4"	6	17.3 lbs
AME2306	13 $\frac{1}{2}$ " x 12 $\frac{1}{2}$ " x 6"	6	20.4 lbs
AM2302-2	Cover—Fits Above Three Trays	6	10.9 lbs
STANDARD WEIGHT (22 GAUGE)			
AME06062	6 $\frac{1}{2}$ " x 6 $\frac{1}{2}$ " x 2 $\frac{1}{2}$ "	6	3.2 lbs
AME06064	6 $\frac{1}{2}$ " x 6 $\frac{1}{2}$ " x 4"	6	5.0 lbs
AME06066	6 $\frac{1}{2}$ " x 6 $\frac{1}{2}$ " x 6"	6	6.6 lbs
HEAVY WEIGHT (20 GAUGE)			
AMS06062	6 $\frac{1}{2}$ " x 6 $\frac{1}{2}$ " x 2 $\frac{1}{2}$ "	6	4.0 lbs
AMS06064	6 $\frac{1}{2}$ " x 6 $\frac{1}{2}$ " x 4"	6	6.0 lbs
AMS06066	6 $\frac{1}{2}$ " x 6 $\frac{1}{2}$ " x 6"	6	7.8 lbs
AM602-2	Cover—Fits Above Six Trays	6	2.4 lbs
STANDARD WEIGHT (22 GAUGE)			
AME10062	10 $\frac{1}{2}$ " x 6 $\frac{1}{2}$ " x 2 $\frac{1}{2}$ "	6	5.3 lbs
AME10064	10 $\frac{1}{2}$ " x 6 $\frac{1}{2}$ " x 4"	6	7.5 lbs
AME10066	10 $\frac{1}{2}$ " x 6 $\frac{1}{2}$ " x 6"	6	9.5 lbs
HEAVY WEIGHT (20 GAUGE)			
AMS10062	10 $\frac{1}{2}$ " x 6 $\frac{1}{2}$ " x 2 $\frac{1}{2}$ "	6	6.0 lbs
AMS10064	10 $\frac{1}{2}$ " x 6 $\frac{1}{2}$ " x 4"	6	8.8 lbs
AMS10066	10 $\frac{1}{2}$ " x 6 $\frac{1}{2}$ " x 6"	6	10.6 lbs
AM402-2	Cover—Fits Above Six Trays	6	3.4 lbs

Order pans and covers separately.



INSTRUMENT TRAYS & FLAT/SOLID COVERS (CON'T)

	DIMENSIONS	CASE	CASE
	L x W x H	LOT	WEIGHT
STANDARD WEIGHT (22 GAUGE)			
AME12062	12 ³ / ₄ " x 6 ⁷ / ₈ " x 2 ¹ / ₂ "	6	6.5 lbs
AME12064	12 ³ / ₄ " x 6 ⁷ / ₈ " x 4"	6	9.3 lbs
AME12066	12 ³ / ₄ " x 6 ⁷ / ₈ " x 6"	6	11.7 lbs
HEAVY WEIGHT (20 GAUGE)			
AMS12062	12 ³ / ₄ " x 6 ⁷ / ₈ " x 2 ¹ / ₂ "	6	7.5 lbs
AMS12064	12 ³ / ₄ " x 6 ⁷ / ₈ " x 4"	6	10.8 lbs
AMS12066	12 ³ / ₄ " x 6 ⁷ / ₈ " x 6"	6	13.5 lbs
AM302-2	Cover—Fits Above Six Trays	6	4.6 lbs
STANDARD WEIGHT (22 GAUGE)			
AME12102	12 ³ / ₄ " x 10 ³ / ₈ " x 2 ¹ / ₂ "	6	8.8 lbs
AME12104	12 ³ / ₄ " x 10 ³ / ₈ " x 4"	6	12.4 lbs
AME12106	12 ³ / ₄ " x 10 ³ / ₈ " x 6"	6	14.8 lbs
HEAVY WEIGHT (20 GAUGE)			
AMS12102	12 ³ / ₄ " x 10 ³ / ₈ " x 2 ¹ / ₂ "	6	10.4 lbs
AMS12104	12 ³ / ₄ " x 10 ³ / ₈ " x 4"	6	14.5 lbs
AMS12106	12 ³ / ₄ " x 10 ³ / ₈ " x 6"	6	17.3 lbs
AM1311-2	Cover—Fits Above Six Trays	6	8.0 lbs
STANDARD WEIGHT (22 GAUGE)			
AME20122	20 ³ / ₄ " x 12 ³ / ₄ " x 2 ¹ / ₂ "	6	17.0 lbs
AME20124	20 ³ / ₄ " x 12 ³ / ₄ " x 4"	6	23.3 lbs
AME20126	20 ³ / ₄ " x 12 ³ / ₄ " x 6"	6	27.7 lbs
HEAVY WEIGHT (20 GAUGE)			
AMS20122	20 ³ / ₄ " x 12 ³ / ₄ " x 2 ¹ / ₂ "	6	20.3 lbs
AMS20124	20 ³ / ₄ " x 12 ³ / ₄ " x 4"	6	26.4 lbs
AMS20126	20 ³ / ₄ " x 12 ³ / ₄ " x 6"	6	33.5 lbs
AM2000-2	Cover—Fits Above Six Trays	6	15.4 lbs

Order pans and covers separately.



PERFORATED INSTRUMENT TRAYS PMP SERIES

- Constructed of heavy duty stainless steel, PMP perforated instrument trays are finished to provide a tough, clean and corrosion-resistant surface for long term use in a variety of applications

	DIMENSIONS	CASE	CASE
	L x W x H	LOT	WEIGHT
PERFORATED MAYO TRAYS			
AM10F-PMT	10" x 6 ¹ / ₂ " x ³ / ₄ "	6	3.5 lbs
AM13F-PMT	13 ¹¹ / ₁₆ " x 9 ¹ / ₁₆ " x ³ / ₄ "	6	9.9 lbs
AM15F-PMT	15 ¹ / ₈ " x 10 ⁵ / ₈ " x ³ / ₄ "	6	10.0 lbs
AM17F-PMT	17" x 11 ⁵ / ₈ " x ³ / ₄ "	6	12.5 lbs
AM19F-PMT	19 ¹ / ₈ " x 12 ³ / ₈ " x ³ / ₄ "	6	14.5 lbs

For excellent steam sterilization without moisture build-up.

	DIMENSIONS	CASE	CASE
	L x W x H	LOT	WEIGHT
HIGH-SIDED PERFORATED TRAYS			
AM952-HPT	8 ¹ / ₈ " x 5" x 2"	12	7.3 lbs
AM1002-HPT	10" x 6 ¹ / ₂ " x 2"	12	10.3 lbs
AM1202-HPT	12 ¹ / ₄ " x 7 ³ / ₄ " x 2 ¹ / ₄ "	12	16.0 lbs
AMS12102-HPT	12 ³ / ₄ " x 10 ³ / ₈ " x 2 ¹ / ₂ "	6	10.4 lbs
AME1650-HPT	16 ¹ / ₂ " x 10" x 2 ¹ / ₂ "	6	12.5 lbs
AMS20122-HPT	20 ³ / ₄ " x 12 ³ / ₄ " x 2 ¹ / ₂ "	6	20.3 lbs

High sides protect fragile instruments during set-up and procedures.

CATHETER TRAYS & DOME COVERS

- Heavy duty and seamless, these trays feature dome covers with strap handles for easy handling

	DIMENSIONS		CASE LOT	CASE WEIGHT
	L	W x H		
AM83	8 $\frac{1}{4}$ "	3 $\frac{1}{4}$ " x 1 $\frac{1}{2}$ "	12	4.5 lbs
AM83-2	Cover—Fits 83 Tray		12	4.0 lbs
AM170-0	17 $\frac{1}{4}$ "	4 $\frac{1}{2}$ " x 2"	12	15.0 lbs
AM170-2	Cover—Fits 170-0 Tray		12	10.3 lbs
AM952	8 $\frac{1}{8}$ "	5" x 2"	12	7.3 lbs
AM953	8 $\frac{1}{8}$ "	5" x 3"	12	8.6 lbs
AM952D-2	Cover—Fits 952 and 953 Trays		12	7.0 lbs

Order pans and covers separately.

STERILIZATION TRAYS/BATH SETS HOT/COLD

- Tray units can be used for either hot or cold sterilization
- Covers feature recessed strap handle for ease and safety in transporting, handling and storing

	DESCRIPTION	DIMENSIONS		CASE LOT	CASE WEIGHT
		L	W x H		
AME20124	Solid Tray	20 $\frac{3}{4}$ "	12 $\frac{3}{4}$ " x 4"	6	23.3 lbs
AMP20122	Perforated Tray	20 $\frac{3}{4}$ "	12 $\frac{3}{4}$ " x 2 $\frac{1}{2}$ "	6	15.4 lbs
AM2000-2	Flat, Solid Cover			6	15.4 lbs
AME20126	Solid Tray	20 $\frac{3}{4}$ "	12 $\frac{3}{4}$ " x 6"	6	27.7 lbs
AMP20124	Perforated Tray	20 $\frac{3}{4}$ "	12 $\frac{3}{4}$ " x 4"	6	14.8 lbs
AM2000-2	Flat, Solid Cover			6	15.4 lbs
AME12104	Solid Tray	12 $\frac{3}{4}$ "	10 $\frac{3}{8}$ " x 4"	6	12.4 lbs
AMP12102	Perforated Tray	12 $\frac{3}{4}$ "	10 $\frac{3}{8}$ " x 2 $\frac{1}{2}$ "	6	8.5 lbs
AM1311-2	Flat, Solid Cover			6	8.0 lbs
AME12106	Solid Tray	12 $\frac{3}{4}$ "	10 $\frac{3}{8}$ " x 6"	6	14.8 lbs
AMP12104	Perforated Tray	12 $\frac{3}{4}$ "	10 $\frac{3}{8}$ " x 4"	6	10.7 lbs
AM1311-2	Flat, Solid Cover			6	8.0 lbs

Order pans and covers separately.

STERILIZATION TRAYS HOT/COLD

- Tray units can be used for either hot or cold sterilization
- Covers feature recessed strap handle for ease and safety in transporting, handling and storing

	CAPACITY	DIMENSIONS		CASE LOT	CASE WEIGHT
		L	W x H		
AME20062	3 $\frac{1}{2}$ qt.	20 $\frac{3}{4}$ "	6 $\frac{3}{8}$ " x 2 $\frac{1}{2}$ "	6	9.9 lbs
AME20064	5 $\frac{1}{4}$ qt.	20 $\frac{3}{4}$ "	6 $\frac{3}{8}$ " x 4"	6	12.4 lbs
AME20066	8 $\frac{3}{4}$ qt.	20 $\frac{3}{4}$ "	6 $\frac{3}{8}$ " x 6"	6	16.5 lbs
AMS20062	3 $\frac{1}{2}$ qt.	20 $\frac{3}{4}$ "	6 $\frac{3}{8}$ " x 2 $\frac{1}{2}$ "	6	11.0 lbs
AMS20064	5 $\frac{1}{4}$ qt.	20 $\frac{3}{4}$ "	6 $\frac{3}{8}$ " x 4"	6	14.6 lbs
AMS20066	8 $\frac{3}{4}$ qt.	20 $\frac{3}{4}$ "	6 $\frac{3}{8}$ " x 6"	6	19.2 lbs
AM2006-2	Flat, Solid Cover—Fits all above pans			6	8.2 lbs

Order pans and covers separately.



MIXING/SOLUTION BOWLS

- Tapered sides and rounded bead design make these basins ideal for storing and pouring

	CAPACITY		DIMENSIONS		CASE	
			DIA. & HT.		LOT	WEIGHT
AM12A	7/8 qt.		6 3/8" x 2 1/2"		12	5.3 lbs
AM18A	1 1/2 qt.		8 3/8" x 3"		12	7.2 lbs
AM22A	3 qt.		9 3/4" x 4"		12	10.8 lbs
AM26A	5 qt.		11 1/8" x 4 1/2"		12	14.0 lbs
AM30A	6 3/4 qt.		12 3/4" x 5 1/4"		6	9.1 lbs
AM32A	8 qt.		13 5/16" x 5 7/16"		6	10.0 lbs
AM40A	13 qt.		16 1/4" x 6"		6	14.0 lbs

SOLUTION BASINS

- Constructed with easy-to-clean large rounded radius corners for easy pouring, and no spin recessed bottom for stability on flat surfaces

	CAPACITY		DIMENSIONS		CASE	
			DIA. & HT.		LOT	WEIGHT
AM130	5 qt.		12" x 4 1/4"		12	18.5 lbs
AM134	7 qt.		13 3/8" x 4 3/8"		6	12.2 lbs
AM136	9 1/2 qt.		14 1/2" x 5"		6	15.1 lbs

EMESIS BASINS

- Vessels feature easy-to-clean seamless design

	CAPACITY		DIMENSIONS		CASE	
			L x W x H		LOT	WEIGHT
AM6	10 oz.		6 5/8" x 3 1/2" x 1 1/2"		12	3.3 lbs
AM8	12 oz.		8 5/8" x 3 3/8" x 1 5/8"		12	4.3 lbs
AM10	26 oz.		9 5/8" x 4 1/2" x 2 1/8"		12	6.9 lbs
AM12	50 oz.		12" x 5 1/2" x 2 3/8"		12	11.0 lbs

IODINE CUPS

	CAPACITY		DIMENSIONS		CASE	
			DIA. & HT.		LOT	WEIGHT
AM6G	6 oz.		3 1/16" x 2"		24	3.5 lbs
AM14G	14 oz.		4 3/8" x 2 5/8"		24	6.8 lbs

SPONGE BOWLS

- Wide mouth ideal for sponges
- Flat bottom and straight sides provide for extra stability in use

	CAPACITY		DIMENSIONS		CASE	
			DIA. & HT.		LOT	WEIGHT
AM52	14 oz.		4 1/4" x 2"		12	3.8 lbs
AM75	20 oz.		5" x 2 3/8"		12	5.1 lbs
AM100	30 oz.		5 7/16" x 2 9/16"		12	5.5 lbs
AM106	38 oz.		6 1/4" x 2 3/4"		12	6.7 lbs
AM148	52 oz.		7 1/8" x 2 3/4"		12	8.0 lbs
AM107	64 oz.		7 9/16" x 3"		12	9.5 lbs



SafeClean™ Instrument and Channel Cleaning Brushes are designed to clean surgical/medical instruments and endoscopic devices. The complete line of SafeClean™ Cleaning Brushes are strong enough to effectively remove surgical residue, yet gentle enough to use on the most delicate instruments and devices. Latex-free.

General Purpose Brushes

General Purpose Brushes can be used throughout all of your instrument and device reprocessing areas to clean uniquely shaped and larger bore surgical, medical, and endoscopic instruments and devices. Latex-free.



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Quantity Per Bag
241007	1mm (0.039")	18"	0.5"	2
241010	1.6mm (0.063")	18"	0.5"	2
241015	1.8mm (0.071")	18"	0.5"	2
241019	2mm (0.079")	12"	1"	2
241020	2.2mm (0.086")	18"	0.5"	2
241024	2.5mm (0.098")	24"	2.5"	2
241025	2.7mm (0.106")	18"	1"	2
241029	3mm (0.118")	12"	1"	2
241030	3.2mm (0.126")	18"	2"	2
241032	4mm (0.158")	12"	1"	2
241033	5mm (0.197")	12"	2.5"	2
241034	5mm (0.197")	24"	2.5"	2
241035	5mm (0.197")	16"	2.5"	2
241040	6mm (0.236")	16"	2.5"	2
241041	6mm (0.236")	24"	2.5"	2
241042	7mm (0.275")	16"	2.5"	2
241043	7mm (0.275")	24"	2.5"	2
241044	10mm (0.393")	16"	2.5"	2
241045	11mm (0.433")	16"	2.5"	2
241046	11mm (0.433")	24"	2.5"	2
241047	10mm (0.393")	24"	2.5"	2
241050	12mm (0.472")	16"	2.5"	2
241055	15mm (0.590")	16"	2.5"	2
241060	20mm (0.787")	16"	2.5"	2

Kits

Make ordering easy with all of the brushes you need in one kit!

Small Cleaning Brush Kit

This kit includes various small general purpose cleaning brushes. Latex-free.

Part No.

Part No.	Quantity
241100	1-Each
241010	
241015	
241020	
241025	
241030	

Large Cleaning Brush Kit

This kit includes various large general purpose cleaning brushes. Latex-free.

Part No.

Part No.	Quantity
241200	1-Each
241035	
241040	
241045	
241050	
241055	
241060	

"Best Sellers" Cleaning Brush Kit

This kit includes a selection of our most popular general purpose cleaning brushes. Latex-free.

Part No.

Part No.	Quantity
241300	1-Each
241010	
241015	
241020	
241025	
241030	
241035	
241040	
241045	
241050	
241055	
241060	

General Cleaning Brush Kit

This kit contains an extensive selection of general purpose brushes, as well as toothbrush and fan tip brushes. Latex-free.

Part No.

Part No.	Quantity
241600	3-Each
241025	
241030	
241035	
241040	
241041	
241045	
241046	
241050	
241055	
241060	
242001	
242002	
242003	
243001	
243004	

Micro/Ophthalmic Brush Kit

This kit contains a selection of our smallest general purpose cleaning brushes. Latex-free.

Part No.

Part No.	Quantity
241700	3-Each
241007	
241010	
241015	
241020	
241025	
241030	

Suction Tube Cleaning Brushes

Suction Tube Cleaning Brushes are designed to fit into and effectively clean the smallest bore devices and instruments used in surgery, MIS, and therapy procedures today. Latex-free.



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Quantity Per Bag
241003	3 French (1mm, 0.039")	6"	0.5"	2
241005	5 French (1.67mm, 0.066")	6"	0.5"	2
241006	6 French (2mm, 0.079")	8"	0.5"	2
241008	7 French (2.3mm, 0.092")	8"	0.5"	2
241009	8 French (2.7mm, 0.105")	8"	0.5"	2
241011	9 French (3mm, 0.118")	8"	0.5"	2
241012	10 French (3.3mm, 0.131")	8"	0.5"	2
241013	11 French (3.7mm, 0.144")	8"	0.5"	2
241014	12 French (4mm, 0.158")	8"	0.5"	2

Fan-Tip Brush


Fan-Tip Brushes feature bristles that wrap around the end of the brush for superior cleaning of bottles, jars, and other large instruments. Latex-free.



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Quantity Per Bag
242001	4.8mm-12.7mm Tapered (0.189"-0.5")	7"	3"	2
242002	13mm (0.512")	6"	3"	2
242003	16mm (0.623")	5"	2.5"	2
242046	9.6mm (0.384")	24"	2.5"	2

Tracheostomy Brushes

Use the right brush for the right job. Our Tracheostomy Brushes are uniquely designed to clean tubes. Latex-free.



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Quantity Per Bag
247001	Fits Trach Sizes 00-SFR 9mm (0.354")	4.75"	1.85"	2
247002	Fits Trach Sizes 6-12FR 13mm (0.512")	8"	3.66"	2

Toothbrush Style

AmeriSurg offers a variety of bristle types, designed to clean all types of instruments and stains including nylon, brass, and steel. Latex-free.



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Quantity Per Bag
243001	13mm (0.512"), Nylon	7.25"	1.5"	2
243002	13mm (0.512"), Brass	7.25"	1.5"	2
243003	13mm (0.512"), Steel	7.25"	1.5"	2

Double-Ended Nylon Cleaning Brushes

Our Double-Ended Cleaning Brush features one end with an extra large brush head for general cleaning and one end with a single row of bristles for more intricate cleaning. Latex-free.



Available
Sterile and
Non-Sterile

Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Quantity
243004	Non-Sterile Short End, 2mm (0.079") Long End, 8mm (0.315")	7"	0.867" 1.38"	2/Bag
243014	Sterile Short End, 2mm (0.079") Long End, 8mm (0.315")	7"	0.867" 1.38"	1/Pack, 12/Packs/Box

Reusable Channel Cleaning Brushes


Channel Brushes with tapered brush heads, elevator sheath, and comfortable round, plastic handles are available in a variety of lengths and materials to fit the range of endoscopic devices. Latex-free.



Part No.	Description	Total Length	Sheath Diameter	Quantity Per Bag
248107	3-5mm	70cm	1.6mm	1
248115	3-5mm	115cm	1.6mm	1
248160	3-5mm	160cm	1.6mm	1
248180	3-5mm	230cm	1.8mm	1
248182	3-5mm — Double Ended	230cm	1.8mm	1
248185	3-5mm — Teflon Coated	230cm	1.8mm	1
248190	3-5mm	120cm	1.8mm	1

Disposable Scope Channel Cleaning Brushes

Disposable Scope Channel Cleaning Brushes have brush heads that range from to 2.5mm – 5mm and will fit scopes as small as 1.1mm. Our disposable scope channel cleaning brushes feature more than double the bristle length than competing low-cost brushes. They have a strong, sturdy handle making it easy to pull through the scope, while still cleaning the scope effectively. Single-use and latex-free.




Disposable Scope Channel Brushes

Part No.	Description	Total Length	Sheath Diameter	Quantity Per Box
248102	2.5mm	230cm	1mm	50
248103	3mm	230cm	1.55mm	50
248104	5mm — Double Ended	230cm	1.85mm	50
248105	5mm	230cm	1.85mm	50

Hand/Nail Brushes

Our Hand/Nail brush is used to clean dirt on hands and nails. It can also be used for general instrument cleaning. Latex-free.



Part No.	Brush Head Diameter	Quantity Per Bag
244001	35mm x 93mm (1.375" x 1.75")	2

Instrument Brush Storage

Convenient storage with our acrylic racks helps organize your cleaning brushes and help ensure you will always have the right brush for the job. Latex-free.

Part No.	Description	Quantity
246001	Brush Rack, Medium (without drip tray)	1
246003	Brush Rack, Large (with drip tray)	1





***NEW Medical and EPA Approved
Antimicrobial Spray***

Tested and proven effective against
a broad range of microbes including:
STAPH, H1N1, E.COLI, MRSA & C-DIFF



**New
product!**

**AMFo2SP
32 oz Spray Bottle**

What is MonoFoil™?

Consumers want antimicrobial products that limit the dangers of potentially harmful bacteria, mold and mildew in their home, work and play environments. MonoFoil™ Antimicrobial Protection is a natural compound that is embedded in or coated on products. MonoFoil™ actually creates a surface that prevents the growth and colonization of microorganisms. This process creates a barrier which inhibits the growth of a broad range of microorganisms including: STAPH, MRSA, E.Coli, Norovirus, and (novel) H1N1, as well as fungi such as “Black Mold.”

Safe, Green and Process Friendly

The MonoFoil™ Technology is a chemically bound, non leaching additive that provides unmatched safety and performance. MonoFoil™ is produced with natural compounds that utilize a organo-functional silane technology that can adhere permanently to any surface or substrate. MonoFoil™ can be applied to surfaces by spraying, dipping, soaking and wiping.


**Ultimate
protection!**



**AMFG128
1 Gallon Container**

How MonoFoil™ Works

MonoFoil's positively charged nitrogen atom ruptures the cell wall of the microbe upon contact. The microbe is then blown apart and completely destroyed, before ever making contact with the surface of the material protected by MonoFoil™.



AMFo2SP 32 Fl oz standard spray container
AMFG128 1 gallon F-Style container
AMF4GBX 4ea gallon box F-Style container

*Effectiveness against an organism should not be interpreted as eliminating, controlling, minimizing or otherwise affecting health conditions which may be associated with the specific organisms.

INSTRUMENT LUBRICANTS

For Surgical Instruments


Product Features:

Helps Prevent Corrosion - Anti-corrosive agents fight corrosion.

Gentle Formula - Formulated with gentle, non-aggressive ingredients.

Extends Instrument Life - Leaves a protective micro-film which prolongs instrument life.

Frees Up Box Locks - Lubricates moving parts of surgical instruments.



AM301

Instrument Lubricant Milk



- AM301-1** One gallon plastic bottle packed four per case
- AM301-5** Five gallon HDPE plastic jug with handle
- AM301-15** Fifteen gallon plastic drum
- AM301-30** Thirty gallon plastic drum
- AM301-55** Fifty-five gallon plastic drum

AM302S

Spray Lubricant



- AM302S** is a convenient, ready to use spray lubricant specially formulated to protect and extend the working life of surgical instruments.
- AM302S-8** 8fl.oz. HDPE sprayer bottles packed 12 per case
- AM302S-16** 16fl.oz. HDPE sprayer bottles packed 8 per case

ENZYMATIC FOAM-SPRAY

For Surgical Instruments

Product Features:

Stops dried on blood - Enzyme foam breaks down & prevents stubborn organic encrustation on surgical instruments.

Eliminates messy spills - No more spillage between the operating room and central cleaning/processing department.

Neutral pH - Safe and effective for stainless steel, tungsten carbide, aluminum and brass instruments.

Bacteriostatic - Microbiological control agents inhibit bacterial growth in the ready-to-use foam solution.

Competitively priced - Costs less without compromising performance.

AM402-22 22 fl.oz. sprayer bottle & sprayer packed twelve per case.

AM402s-Gal One gallon bottle packed four per case. (refill for 22fl.oz. bottle)



AM402
Enzymatic Foam Spray

MEDI-SHEEN

Surgical Instrument Reconditioner

Product Features:

Tough on stains, gentle on instruments - Medi-Sheen is a powerful stain and rust remover coupled with corrosion inhibitors to recondition and provide lasting protection to stainless steel instruments.

Frees up frozen box locks - In most cases, instruments with frozen box locks that are out of service are restored with Medi-Sheen.

Ideal for cleaning steam autoclaves - Medi-Sheen's robust stain and rust removing agents quickly cut through scale and deposit, typical of steam sterilizers.

Cost Competitive - Medi-Sheen gets the job done faster for less, compared to similar competitive products.

AM800-1 One gallon plastic bottle packed four per case

AM800-5 Five gallon plastic jug with handle

AM800-15 Fifteen gallon plastic drum

AM800-30 Thirty gallon plastic drum



AM800
Medi-Sheen

Alton Dean Pressure Infusor (Bag Squeezer)

When connected to a pneumatic pressure source or to an electrical outlet, Alton Dean Pressure Infusors automatically maintain pressure on blood and IV fluid bags. Pressure range is zero to 300mm HG. An optional Quick Disconnect Coupler allows easy separation of the pressure hose from the Infusor. An optional Hand Pump maintains pressure during transport or when away from a pressure source.

When you get tired of hand pumping...

We're ready to start pumping for you.



ALTON DEAN AUTOMATIC PRESSURE INFUSOR NOW AVAILABLE IN ELECTRIC & PNEUMATIC

- **AUTOMATIC - NO DROP IN PRESSURE**
Maintains any selected pressure from 0 to 300 millimeters of Mercury (0-5.8 psi) on a 500 ml bag.
- **NO HAND PUMPING** - Plugs into oxygen or air outlets for automatic pressure maintenance.
- **FAST LOADING** - Hinged, transparent door allows fast loading of fluid bag and unobstructed view of fluid level.
- **LIGHTWEIGHT** - Durable materials.
- **COMPACT** - Hangs on any I.V. pole.
- **RELIABLE** - Non-electric components.
- **COST EFFECTIVE** - Built for long life, and backed with a one-year warranty.
- **UNIVERSAL** - For the pressure infusion needs of every department in the hospital.

Pressure Infusor
Product #70-5003 Pneumatic
Product #70-5013 Electric

Trump-It Suction Irrigation Set

- No Batteries
- Dual Spike
- Tubing Split (Irrigation/Suction)
- Multiple Probes & Cautery Probes
- Ergonomic Hand Piece
- Snap on Probes



70-5001
Suction Irrigation Set

LATEX FREE & DEHP FREE

Lapro-Flo Insufflation Tubing

with .1 micron ULPA Filter 10" of DEHP FREE Soft Tubing

Insufflation Filter Sets that allow for high volume at low pressure settings.

Filters Particulate Matter of .1 Micron With 99.997% Efficiency

Hydrophobic ULPA Filtration

**"SNAP-ON"
Connector**



70-5006
O-Ring Connector
Smith & Nephew



70-5007
15/22mm Connector
Wolf, Storz, Linvatec



70-5008
Funnel Connector
Stryker



70-5009
Luer Lock Connector
Dual Luer

Insufflation Tubing Set

works with Storz Heated Insufflator

This Insufflation Tubing Set features:

- 1/2" Tubing
- 0.01 Micron Filter, 99.999% Efficient
- Hydrophobic ULPA Filter
- DEHP Free
- 11' of Soft High Flow Tubing



70-5010
"Heat Activated"
Insufflation Tubing Set

Amerisurg Trump-It Accessories



- 70-4000** Spatula Probe, 33cm
- 70-4001** L-Hook Probe, 33cm
- 70-4002** J-Hook Probe, 33cm
- 70-4003** Poole Suction, 10mm/33cm
- 70-4004** Poole Suction, 5mm/33cm
- 70-4005** Probe w/Holes, 10mm/33cm
- 70-4006** Probe w/Holes, 5mm/33cm
- 70-4007** Probe w/Holes, 5mm/45cm
- 70-4008** Micro Dissector 5mm/2.5/33cm



Laparoscopic Kittner
Blunt Dissecting End 20/case
Product #70-5004



Laparoscopic Starter Kit
1ea Insufflation Tubing with Filter
1ea Anti-Fog Solution 6mm Bottle & Sponge
Product #70-5005

**Eliminate
Scope Condensation**

Scope Warmer
Product #70-5000
Box of 10/Case of 30

SCOPE WARMER



Our new longer lasting, double sided scopewarmer pre-warms endoscopes to body temperature to keep them free from condensation. It's easy to use and each side can provide over one hour of warming time. No batteries or difficult to activate charcoal-type warmers. No need for costly decontamination, sterilization, and additional sterile water for thermos-type warmers.

- Disposable, sterile heating device
- Stays warm for over 2 hours
- Warms endoscopes in less than a minute
- Inexpensive and easy to use
- Latex and DEHP Free
- Environmentally friendly
- Contains Saline only - no chemicals

**NEW & IMPROVED
EXTENDED LIFE
Warms Up
To 2 Hours**




Anti-Fog Solution
6mm Sterile Bottle & Sponge 20/case
Product #70-5011



Verres Insufflation Needle
12 Gauge, 12mm
Product #70-5002

Endoscopic Fiberoptic Light Cables



10 Colors Available:

Gray	
Green	
Yellow	
Opaque	
Red	
Orange	
White	
Blue	
Pink	
Purple	












How to Order Your Endoscopic Cable












1. Select your light source end
2. Select your instrument end fitting
3. Select your cable color

Light Source End


Instrument End

.55 NUMERICAL APERTURE

A		ACMI
AO		ACMI/OLYMPUS
S		STORZ
W		WOLF
WX		EXTENDED WOLF
U		UNIVERSAL *
P		PILLING
PA		PENTAX
O		OLYMPUS
M		MACHIDA
D		DESIGN FOR VISION

ACMI MALE		A
ACMI FEMALE SNAP ON		F
WOLF FEMALE SNAP ON		W
*UNIVERSAL INSTRUMENT		U
EXTERNAL WOLF CLIP		EXW
PILLING BAYONET		P
STORZ/OLYMPUS		S
OLYMPUS RETRACTABLE		ORF
ACMI FEMALE CLIP		FC
DESIGN FOR VISION		D
ACMI FEMALE SWIVEL		FS

Urology Resecting Cutting Loops



ORDERING INFORMATION

Olympus ACMI Classic Resectoscope

Order Number	Resection / Ablation Device	Units/ Box
ACC-CL2412	Cutting Loop 24 Fr. (.012 wire)	12
ACC-CL2612	Cutting Loop 26 Fr. (.012 wire)	12
ACC-CL2812	Cutting Loop 28 Fr. (.012 wire)	12
ACC-CL2415	Cutting Loop 24 Fr. (.015 wire)	12
ACC-CL2615	Cutting Loop 26 Fr. (.015 wire)	12
ACC-CL2815	Cutting Loop 28 Fr. (.015 wire)	12
ACC-CP2400	OMNI-LOOP Resector 24 Fr.	6
ACC-CP2600	OMNI-LOOP Resector 26 Fr.	6
ACC-CP2800	OMNI-LOOP Resector 28 Fr.	6
ACC-RB243M	McCausland 3mm Roller-ball	6
ACC-RB245M	5mm Roller-ball	6
ACC-RL243M	3mm Roller Barrel	6
ACC-DV243M	3mm Dimpled Vaporization Roller	6
ACC-CK2400	Collings Knife 24 Fr.	6

Olympus ACMI USA/Elite Resectoscope

Order Number	Resection / Ablation Device	Units/ Box
ACU-CL2412	Cutting Loop 24 Fr. (.012 wire)	12
ACU-CL2612	Cutting Loop 26 Fr. (.012 wire)	12
ACU-CL2812	Cutting Loop 28 Fr. (.012 wire)	12
ACU-CL2415	Cutting Loop 24 Fr. (.015 wire)	12
ACU-CL2615	Cutting Loop 26 Fr. (.015 wire)	12
ACU-CL2815	Cutting Loop 28 Fr. (.015 wire)	12
ACU-CP2400	OMNI-LOOP Resector 24 Fr.	6
ACU-CP2600	OMNI-LOOP Resector 26 Fr.	6
ACU-CP2800	OMNI-LOOP Resector 28 Fr.	6
ACU-RB243M	McCausland 3mm Roller-ball	6
ACU-RB245M	5mm Roller-ball	6
ACU-RL243M	3mm Roller Barrel	6
ACU-DV243M	3mm Dimpled Vaporization Roller	6
ACU-CK2400	Collings Knife 24 Fr.	6

Storz Single Stem Resectoscope

Order Number	Resection / Ablation Device	Units/ Box
KS1-CL2412	Cutting Loop 24 Fr. (.012 wire)	6
KS1-CL2712	Cutting Loop 27 Fr. (.012 wire)	6
KS1-CL2414	Cutting Loop 24 Fr. (.014 wire)	6
KS1-CL2714	Cutting Loop 27 Fr. (.014 wire)	6
KS1-CP2400	OMNI-LOOP Resector 24 Fr.	6
KS1-CP2700	OMNI-LOOP Resector 27 Fr.	6
KS1-RB243M	McCausland 3mm Roller-ball	6
KS1-RB245M	5mm Roller-ball	6
KS1-RL243M	3mm Roller Barrel	6
KS1-DV243M	3mm Dimpled Vaporization Roller	6
KS1-CK2400	Collings Knife 24 Fr.	6

Storz Double Stem Resectoscope

Order Number	Resection / Ablation Device	Units/ Box
KS2-CL2412	Cutting Loop 24 Fr. (.012 wire)	6
KS2-CL2712	Cutting Loop 27 Fr. (.012 wire)	6
KS2-CL2414	Cutting Loop 24 Fr. (.014 wire)	6
KS2-CL2714	Cutting Loop 27 Fr. (.014 wire)	6
KS2-CP2400	OMNI-LOOP Resector 24 Fr.	6
KS2-CP2700	OMNI-LOOP Resector 27 Fr.	6
KS2-RB243M	McCausland 3mm Roller-ball	6
KS2-RB245M	5mm Roller-ball	6
KS2-RL243M	3mm Roller Barrel	6
KS2-DV243M	3mm Dimpled Vaporization Roller	6
KS2-CK2400	Collings Knife 24 Fr.	6

Olympus OES Resectoscope

Order Number	Resection / Ablation Device	Units/ Box
OES-CL2412	Cutting Loop 24 Fr. (.012 wire)	12
OES-CL2612	Cutting Loop 26 Fr. (.012 wire)	12
OES-CL2812	Cutting Loop 28 Fr. (.012 wire)	12
OES-CL2414	Cutting Loop 24 Fr. (.014 wire)	12
OES-CL2614	Cutting Loop 26 Fr. (.014 wire)	12
OES-CL2814	Cutting Loop 28 Fr. (.014 wire)	12
OES-CP2400	OMNI-LOOP Resector 24 Fr.	6
OES-CP2600	OMNI-LOOP Resector 26 Fr.	6
OES-CP2800	OMNI-LOOP Resector 28 Fr.	6
OES-RB243M	McCausland 3mm Roller-ball	6
OES-RB245M	5mm Roller-ball	6
OES-RL243M	3mm Roller Barrel	6
OES-DV243M	3mm Dimpled Vaporization Roller	6
OES-CK2400	Collings Knife 24 Fr.	6



OMNI-LOOP
Coagulating Resector

Each Unit is in a protective blister inside a Tyvek to film pouch, STERILE, packaged in Dispenser Box.



Made in U.S.A.

LaserGuard™ Technology

This patented accessory for laser fibers:

- Indicates fiber failure between tip of fiber and connector to laser
- Prevents scope damage due to malfunction of laser fiber
- Notifies user of risk of damage to accessories used in stone removal

ScopeSafe™ Technology

- Filters errant energy that causes scope damage



LaserGuard™ and ScopeSafe™ are registered trademarks of Optical Integrity, Inc. and are protected by patents: 7,540,668 7,204,645 6,948,862 and 7,090,411 and patent pendings.



Let's Define the Problem

Endoscopes are damaged by laser fibers one of three ways:

1. Fiber breaks inside the scope while lasing
2. Fiber pulled into the scope while lasing
3. Mechanical damage from fiber insertion

Optical Integrity Answer to #1: ScopeSafe™

Fibers break while lasing because there is too much energy in the cladding. Optical Integrity Inc. designed, manufactured and patented new optical filtering technology that prevents laser energy from getting into the cladding in the first place.




Optical Integrity Answer to #2: LaserGuard™

Optical Integrity designed, manufactured, and patented the LaserGuard™ to shut down the laser when the fiber is pulled back into the scope while lasing. Also prevents laser damage to baskets and grabbers in the surgical field of view.

Optical Integrity Answer to #3: LaserGuard Sheath™

Optical Integrity's patented LaserGuard Sheath™ prevents mechanical damage to endoscopes from fiber insertion.

1. Does not affect the deflection of scope.
2. Does not reduce water flow.
3. Allows fiber insertion into deflected scope.



Comparison of Optical Integrity ScopeSafe™ Fibers with LaserGuard™ to Competitor's Fibers

	OI	Competitors	
Smallest Size (µm core)	145	150	272
Blast Shield Warranty	Yes	No	No
Scope Warranty W/LaserGuard™	Yes	No	No
Scope Pull-in Protection regardless of field clarity or laser settings	Yes	No	No
Scope Fiber Break Protection	Yes	No	No
Basket Protection, Stone Cone or Guide wire laser protection	Yes	No	No
Tip Degradation for both side fire and flat tip fibers	Yes	No	No
Allows for fiber insertion into a fully deflected scope	Yes	No	No
Detect end life of fibers	Yes	No	No



Model Number	Description
--------------	-------------

ScopeSafe™ Laser Fiber with Stone Grasper

Sterile, Single-use (**Disposable**):

10547	145m Core Polyamide Fiber (for use with Stone Grasper below)
10558	2.4 French 115cm Stone Grasper w/ Laser Fiber Channel (Fiber # 10547)

ScopeSafe™ Laser Fibers for Holmium Lasers with SMA 905 Port

Sterile, Single-use (**Disposable**) Fibers:

P0012-0200	200µm Core ScopeSafe™ Fiber, Single-use
P0012-0272	272µm Core ScopeSafe™ Fiber, Single-use
P0012-0300	300µm Core ScopeSafe™ Fiber, Single-use
P0012-0365	365µm Core ScopeSafe™ Fiber, Single-use
P0012-0550	550µm Core ScopeSafe™ Fiber, Single-use
P0012-1000	1000µm Core ScopeSafe™ Fiber, Single-use

LaserGuard™ Technology for Holmium Lasers with SMA 905 Port

Sterile, Single-use Fibers (**Disposable**), to be used only with LaserGuard™ technology.

10533	200µm LaserGuard™ Fiber, Single-use
10534	272µm LaserGuard™ Fiber, Single-use
10535	320µm LaserGuard™ Fiber, Single-use

LaserGuard™ Device (Accessory for use with LaserGuard™ Fibers):

1047x	LaserGuard™ Device
-------	--------------------

Fiber Accessories

Fiber Stripping Tools:

10305	Stripper Tool for 200µm Fiber
10306	Stripper Tool for 272µm Fiber
10307	Stripper Tool for 300µm Fiber
10557	Stripper Tool for 320µm Fiber
10308	Stripper Tool for 365µm Fiber
10309	Stripper Tool for 550µm Fiber
10310	Stripper Tool for 1000µm Fiber

Fiber Cutting and Cleaving Tools:

10554	Ceramic Scissors
10095	Ceramic Cleaver

Laser Eye Ware:

10555	Laser Glasses for 2100nm (Holmium) Wavelength
10556	Laser Goggles for 2100nm (Holmium) Wavelength

Custom Fiber Assemblies

We can also custom design/build products to your specific needs and specific laser systems.

OB/GYN Catheters/Disposables



AMI1151 Zinnanti Uterine Manipulator/Injector

- Pear shaped balloon for optimal hold which is a better seal at the internal os
- Reduced catheter shaft (4.5mm) allows for easy insertion
- Flexible extension prevents cracking of luer fitting

LATEX FREE



AMI1171 EZ Glide Manipulator

- New sliding handle adjusts after insertion and then locks into place to maintain a secure trauma free hold on the uterus
- Handle has larger grip for better leverage during manipulation
- Pear shaped balloon provides optimal hold which forms a better seal at the internal os
- Flexible extension prevents cracking of luer fitting
- Longer shaft to better accommodate Plus-Size patients
- Central lumen provides for easy injection
- Atraumatic beveled catheter tip protects uterine wall


LATEX FREE

Ideal for:

- Laparoscopic Tubal Occlusion
- Minilaparotomy
- Fertility Exams
- Laparoscopic Assisted Vaginal Hysterectomy (LAVH)
- Diagnostic and Operative Laparoscopy
- Salpingoplasty
- Laparoscopic Supracervical Hysterectomy (LSH)

Uterine Dilators

Reusable Teflon Dilators - Autoclavable



BENEFITS:

- Autoclavable
- Gas sterilizable
- Soakable
- Less traumatic to patient with stenotic cervix

FEATURES:
























- Made of Teflon with longevity of stainless
- Self-lubricating, malleable material
- Three different tip lengths

Dilator set is used sequentially in the order of rings around the handle.

OS 1
CANAL 2
FUNDUS 3

- AM 01-100** Cervical Os Finder, Single Ended, Diameter 1.75mm to 3.8mm, Depth 2.3cm, Length 20cm (8")
- AM 01-200** Uterine Canal Finder, Single Ended, Diameter 1.75mm to 3.8mm, Depth 4.5cm, Length 20cm (8")
- AM 01-300** Uterine Cavity Finder, Single Ended, Diameter 1.75mm to 3.8mm, Depth 8.5cm, Length 20cm (8")
- AM 01-500** Flexible Dilator Set of 3, Single Ended, Diameter 1.75mm to 3.8mm, Depth 2.3cm to 8.5cm, Length 20cm (8")

Loop Electrodes

	Ball	AM13-003	3mm, Shaft: 13cm (5")
	Square	AM14-005	5mm, Shaft: 13cm (5")
		AM29-505	5mm x 5mm, Square, Shaft: 13cm (5")
		AM06-101	10mm x 10mm, Square, Shaft: 13cm (5")
	Loop	AM06-108	10mm x 8mm, Square, Shaft: 13cm (5")
		AM03-505	5mm x 5mm, Shaft: 13cm (5")
		AM04-105	10mm x 5mm, Shaft: 13cm (5")
		AM22-107	10mm x 7mm, Shaft: 13cm (5")
		AM28-108	10mm x 8mm, Shaft: 13cm (5")
		AM05-101	10mm x 10mm, Shaft: 13cm (5")
		AM23-155	15mm x 5mm, Shaft: 13cm (5")
		AM27-157	15mm x 7mm, Shaft: 13cm (5")
		AM07-158	15mm x 8mm, Shaft: 13cm (5")
		AM08-151	15mm x 15mm, Shaft: 13cm (5")
		AM08-153	15mm x 13mm, Shaft: 13cm (5")
		AM09-208	20mm x 8mm, Shaft: 13cm (5")
		AM24-201	20mm x 10mm, Shaft: 13cm (5")
		AM28-212	20mm x 12mm, Shaft: 13cm (5")
		AM10-203	20mm x 13mm, Shaft: 13cm (5")
		AM11-205	20mm x 15mm, Shaft: 13cm (5")
		AM25-202	20mm x 20mm, Shaft: 13cm (5")
		AM12-251	25mm x 10mm, Shaft: 13cm (5")
		AM26-150	15mm x 10mm, Shaft: 13cm (5")

Sterile, Disposable, Each Sold as a Box of 5

- AM02-012** Conductive Adhesive Electrode Pads, Disposable, Adult; Box of 12
- AM03-025** Conductive Adhesive Electrode Pads, Disposable for ZSI ESU; Box of 25
- AM04-025** Conductive Adhesive Electrode Pads, Dual Dispersive Disposable; Box of 25
- AM39-025** Conductive Electrode Pads With Cable, Disposable, Cooper; Box of 25
- AM40-050** Conductive Electrode Pads Without Cable, Disposable; Box of 50
For Cooper 1000 and 6000 ESU
- AM05-020** Hand Switching Pencil and Cord, Disposable; Box of 20
- AM09-020** Foot Control Handle and Cord #950, Disposable; Box of 20 For Cooper
- AM20-300** Smoke Evacuator Filters for MiniVac; Case of 3
- AM21-375** Smoke Evacuation Tubing Assembly; Corrugated Tubing, Reducer,
3/8" Tubing, Box of 3, Assemblies
- AM22-006** Smoke Evacuator Tubing, 3/8" Diameter, Case of 6, 61cm (24")
- AM23-100** Smoke Evacuator Accessory Kit, 1 Filter and 1 Tubing Assembly
- AM26-500** DSE Tubing, Box of 50
- AM31-012** Monsel's Solution, 8ml, with Applicator, Package of 12
- AM32-012** Lugol's Iodine Solution, 5ml Doses, Package of 12
- AM38-101** Versa Vac Smoke Evacuator Filter
- AM34-078** Corrigated Tubing, 7/8" x 6'
- AM36-010** Reducer for Smoke Evacuator, Disposable, Single Patient Use,
Connects Pre-Filter to Patient Tubing Set, 3/8", Box of 10
- AM36-100** Pre-Filter for Smoke Evacuator Replacement Filter for Cooper System 6080, Box of 10
- AM35-001** Reducer for Smoke Evacuator, Connects Pre-Filter (AM-36-100) to
Corrigated Tubing and "s" Tubing
- AM33-010** Patient Tubing Set, 6' Disposable Tubing with DSE Tube, 10 Sets per Box



AM05-020
Hand Switching Pencil

AM23-100
Smoke Evacuator Accessory Kit

AM02-012
Conductive Adhesive Electrode Pads

Disposable Electrodes

All disposable electrodes utilize the standard 3/32" (2.3 mm) stainless steel shaft. All molded electrodes feature safety grip insulators combining patient and user safety with easy insertion into and removal from the surgical pencil. All disposable electrodes are manufactured to the highest standards. They come individually packaged sterile.



Dermal Tips

AMA804 – Sharp dermal tip – non-sterile (not packaged individually). 100/box



AMA805 – Sharp dermal tip – sterile (packaged individually). 50/box



AMA806 – Blunt dermal tip – non-sterile (not packaged individually). 100/box



AMA807 – Blunt dermal tip – sterile (packaged individually). 50/box



Blades

AMES01 standard blade - 50/box



AMES04 extended blade - 25/box



AMES37 modified blade - 25/box



AMES39 extended modified blade - 25/box



AMES18 angled blade - 5/box



AMES54 extended modified blade 4" - 25/box



AMES55 extended blade 4" - 25/box



Balls

AMES06 extended 3 mm ball - 5/box



AMES21 5/32" ball - 25/box



Needles

AMES02 standard needle - 25/box



AMES03 extended needle - 25/box



AMES38 modified needle - 25/box



AMES40 extended modified needle - 25/box



Bovie Cauteries

Low-Temperature Cauteries

The Best Cauteries Ever Made

Aaron manufactures four different cauteries specifically designed for ophthalmology. These patented Aaron RAM™ cauteries have an impressive four-year shelf-life and are quite simply the best cauteries ever made. With a variety of low temperatures and fine tips the surgeon is assured of precise pinpoint hemostasis. All disposable cauteries are individually packaged sterile 10 per box.


AMAA00
low-temp fine tip
1300° F (704° C) - 10/box




AMAA02
low-temp adjustable fine tip
700°-1200° F (371°-649° C) - 10/box




AMAA04
low-temp elongated fine tip
1100° F (593° C) - 10/box




AMAA90
low-temp micro fine tip
850° F (454° C) - 10/box
**AMAA90 has one-year shelf-life due to AAA batteries.*



High-Temperature Cauteries

Subungual Hematoma

Cautery

The standard procedure for relief of subungual hematoma is the use of a fine tip high-temp cautery. This procedure gives virtually instant relief.

AMAA01
high-temp fine tip 2200° F (1204° C) - 10/box



AMAA01

Nail Drill Bit

Also available is the Aaron Bit, replacement drill bit.

AM0007
Aaron bit, replacement drill tip - 10/box



AM0007

High-Temperature Cauteries

The Best Cauteries Ever Made

Aaron manufactures more cauteries than anyone in the world. Therefore, Aaron has the ability to put the right cautery in the physician's hand regardless of the procedure. High temperature cauteries are popular for many types of procedures including, pinpoint hemostasis, evacuation of a subungual hematoma, controlling diffuse bleeders, vasectomies, and sculpting woven grafts. Individually packaged sterile, 10 per box.


AMAA01

high-temp fine tip 2200° F (1204° C) - 10/box




AMAA03

high-temp loop tip 2200° F (1204° C) - 10/box




AMAA05

high-temp loop tip w/extended 2" shaft
2200° F (1204° C) - 10/box




AMAA09

high-temp loop tip w/extended 5" shaft
2200° F (1204° C) - 10/box




AMAA11

high-temp adjustable fine tip
1500-2100° F (816°-1149° C) - 10/box




AMAA17

high-temp fine tip w/extended 2" shaft
2200° F (1204° C) - 10/box




AMAA21

high-temp vasectomy tip
1800° F (982° C) - 10/box




AMAA25

high-temp fine tip micro
2000° F (1093° C) - 10/box




AMAA27

high-temp loop tip micro
2000° F (1093° C) - 10/box




AMAA29

high-temp vasectomy tip micro
1600° F (871° C) - 10/box



Ophthalmic Burr Handle

AM0010

Aaron Burr 2 includes on disposable 1 mm burr tip [AM0001] for corneal rust ring removal, "AA" alkaline battery and a foam-lined case, a diamond burr tip [AM0011] may be purchased separately for use in polishing the pterygium bed after surgical removal and for lid margin lesions - 1/box

AM0010 - Aaron Burr 2



Ophthalmic Burr Tips

AM0011

diamond burr 2.5 mm for use in conjunction with the AM0010 or equivalent - 1/box

AM0012

diamond burr 5 mm for use in conjunction with the AM0010 or equivalent - 1/box

AM0001

1 mm burr tip for use with the AM0010 or equivalent for corneal rust ring removal - 10/box

Burr Tips



AMAB01

1 mm burr tip for use with AMAB03, AMAB04 or equivalent - 10/box

AMAB05

1/2 mm burr tip for use with AMAB03, AMAB04 or equivalent - 10/box

Neuro-Pulse™ Description

Neuro-Pulse™ nerve locators were designed for identifying exposed motor nerves while at the same time reducing the possibility of accidental nerve damage or severence. These self powered units are portable and battery operated. Ideal for head, neck, hand, plastic, and facial applications.

AM0003 (Red LED)

Neuro-Pulse™ has three output settings:
1/2 mA, 1 mA, 2 mA at 3 volts DC - 10/box

AM0003Y (Yellow LED, for use in Europe)

Neuro-Pulse™ has three output settings:
1/2 mA, 1 mA, 2 mA at 3 volts DC - 10/box



Ultracell® Eye Products

American Surgical distributes a complete line of sponge products used for management of fluids during ophthalmic procedures. Our PVA and cellulose sponge products are used during cataract and other refractive surgeries. For LASIK procedures, our UltraSpear™, ultra clean PVA surgical eye spears, are ideal for smoothing the flap and manipulating tissue. We also offer corneal light shields, drains, and wicks in various sizes and capacities.

Specially Formulated PVA Material:


- Ultracell® PVA is a synthetic sponge composed of polyvinyl alcohol
- Specially formulated PVA material is virtually lint-free
- Made of ultra fine pores
- Continuous open pore construction creates strength as well as softness
- Highly absorbent
- Instantaneous expansion
- Biocompatible

LASIK Products

American Surgical offers products specially formulated for use during LASIK surgery. We recommend the following line of products for use during LASIK surgery.

LASIK PVA Spears and Disk

American Surgical offers UltraSpear™ LASIK eye spears. The UltraSpear™ material features our exclusive patent-pending particulate extraction process to make it the cleanest eye spear we offer. This means they're virtually particulate free. UltraSpear™ Classic is a traditional compressed shape and is the fastest expanding spear available. Our pre-expanded UltraSpear™ features instantaneous wicking. The pre-expanded tip allows surgeons to adjust the firmness for better control during tissue manipulation. Our LASIK Disk is used for blotting and drying the corneal flap and bed. Latex-free and intended for single-use.



Part No.	Description	Quantity
AM40800	UltraSpear™ Expanded	2/Pack, 20/Box
AM40805	UltraSpear™ Expanded	10/Pack, 10/Box
AM40813	UltraSpear™ Classic	3/Pack, 20/Box
AM40815	UltraSpear™ Classic	10/Pack, 10/Box
AM40804	LASIK Disk	1/Pack, 20/Box

LASIK Drains & Shields

American Surgical also offers a family of other Ultracell® brand PVA products specially designed for the fluid management needs unique to LASIK surgery. Products include drains and shields for both nasal and superior hinges. Latex-free and intended for single-use.



Part No.	Description	Quantity
AM40820	LASIK Shield (7mm flat spot)	1/Pack, 20/Box
AM40821	LASIK Shield (6mm flat spot)	1/Pack, 20/Box
AM40822	LASIK Shield (4mm & 8mm flat spots)	1/Pack, 20/Box
AM40830	LASIK Drain (11.5mm opening)	1/Pack, 20/Box
AM40831	LASIK Drain (10.5mm opening)	1/Pack, 20/Box

Corneal Marking Pad

American Surgical's Corneal Marking Pad can be used with many refractive/corneal marking instruments. The pre-inked, ready-to-use gentian violet marking pad saves prep time and allows you to focus on your patient. Sterile, single-use, and latex-free.

Part No.	Description	Quantity
AM40490	Corneal Marking Pad	1/Pack, 10/Box




Cataract and Refractive Surgery Products

American Surgical provides a complete line of Ultracell® brand surgical eye spears including both cellulose and PVA spears used in the rapid recovery of fluids during cataract and other refractive surgeries.

PVA & Cellulose Eye Spears

American Surgical distributes both PVA and cellulose eye spears for use in ophthalmic surgery. The cellulose spear is the standard spear used in cataract surgery. The PVA eye spear is used in both cataract and LASIK surgeries and has less particulate than the cellulose spear. Latex-free and intended for single-use.



Part No.	Description	Quantity
AM40400-8	PVA Spears	5/Pack, 20/Box
AM40405-8	PVA Spears	10/Pack, 20/Box
AM40410	Cellulose Spears	5/Pack, 20/Box
AM40413	Cellulose Spears	6/Pack, 24/Box
AM40415	Cellulose Spears	10/Pack, 20/Box

Part No.	Description	Quantity
AM40401	PVA Points	5/Pack, 20/Box
AM40406	PVA Points	10/Pack, 20/Box
AM40411	Cellulose Points	5/Pack, 20/Box
AM40416	Cellulose Points	10/Pack, 20/Box
AM40412	Cellulose Points and Strips	2/Pack, 6/Pack, 20/Box

Other Cataract and Refractive Surgery Products

American Surgical provides a variety of Ultracell® PVA products including eye drains, wicks, shields, ring, and pledgets to aid in the continuous collection and control of excessive fluids during surgery. Drains include a fluid collection bag, while wicks are the PVA strip only. Latex-free and intended for single-use.

Part No.	Description	Quantity
AM40420	Corneal Light Shield (7mm)	1/Pack, 20/Box
AM40421	Corneal Light Shield (8mm)	1/Pack, 20/Box
AM40422	Corneal Light Shield (9mm)	1/Pack, 20/Box
AM40430	Eye Drain (80cc capacity)	1/Pack, 20/Box
AM40431	Eye Wick (4mm x 170mm)	2/Pack, 10/Box
AM40435	Eye Drain (400cc capacity)	1/Pack, 10/Box
AM40440	Bloomberg Ophthalmic Ring	1/Pack, 20/Box
AM40450	Eye Kit (1 Wipe/Wick, 1 Light Shield)	1/Pack, 20/Box
AM40451	Cataract Kit (5C. Spear, Lgt Shld, Eye Dm, Inrt Wp)	1/Pack, 20/Box
AM40825	Eye Pledget	10/Pack, 20/Box



Ultracell® ENT Products

American Surgical distributes a complete line of medical sponge products for otolaryngology including the nasal, sinus, epistaxis, and otology sub-specialties. Along with our products for instrument cleaning and care, all Ultracell® ENT sponge products are produced using specially formulated PVA material making them both highly absorbent and soft. For the surgeon, this means high performance. For the patient, this means increased comfort.

Nasal Packs

American Surgical provides a complete line of Ultracell® brand PVA sponges for postoperative packing following nasal surgery.

- Used following septoplasty, rhinoplasty, and turbinectomy procedures
- May aid in holding septum firmly in place
- Available in a range of sizes: 3.0cm pediatric to 10.0cm
- Select sizes available in an anatomical shape
- Popular sizes also available in our patented layered design
- Excellent alternative to cumbersome gauze packing
- Strung for safety
- Latex-free
- Single-use



Regular Nasal

Part No.	Description	Quantity
AM10096-C	3cm Pediatric Nasal Slim Pack	1/Pack, 20/Box
AM10100-C	4.5cm Nasal Pack	1/Pack, 20/Box
AM10101-C	4.5cm Nasal Slim Pack	1/Pack, 20/Box
AM10105-C	6.0cm Nasal Pack	1/Pack, 20/Box
AM10106-C	6.0cm Nasal Slim Pack	1/Pack, 20/Box
AM10110-C	8.0cm Nasal Pack	1/Pack, 20/Box
AM10111-C	8.0cm Nasal Slim Pack	1/Pack, 20/Box
AM10115-C	10.0cm Nasal Pack	1/Pack, 20/Box
AM10120-C	8.0cm Anatomical Pack	1/Pack, 20/Box
AM10125-C	8.0cm Anatomical Slim Pack	1/Pack, 20/Box
AM10126-C	5.5cm Pediatric Anatomical Pack	1/Pack, 20/Box
AM10127-C	5.5cm Pediatric Anatomical Slim Pack	1/Pack, 20/Box



Layered Nasal

Part No.	Description	Quantity
AM10100-L	4.5cm Nasal Pack, Layered	1/Pack, 20/Box
AM10105-L	6.0cm Nasal Pack, Layered	1/Pack, 20/Box
AM10110-L	8.0cm Nasal Pack, Layered	1/Pack, 20/Box
AM10115-L	10.0cm Nasal Pack, Layered	1/Pack, 20/Box
AM10120-L	8.0cm Anatomical Pack Layered	1/Pack, 20/Box

UltraPlus 5000™ Nasal Pack

American UltraPlus 5000™ Nasal Pack features an 8.0cm nasal pack inside custom wrapping to help ease the pain of removing the nasal pack. Latex-free, intended for single-use, and strung for safety.



Part No.	Description	Quantity
AM50110	UltraPlus 5000 8.0cm Nasal Pack	1/Pack, 10/Box
AM50111	UltraPlus 5000 8.0cm Nasal Slim Pack	1/Pack, 10/Box

Sinus Packs

American Surgical offers a complete line of Ultracell® brand sinus packing specifically designed for use following sinus surgeries.

- Ideal packing following FESS (Functional Endoscopic Sinus Surgery)
- Used for post-operative packing or dressing of sinus cavity
- Available in a range of lengths and thicknesses
- Patented layered design allows for easy manipulation
- Strung for safety
- Latex-free
- Single-use



Part No.	Description	Quantity
AM20200-C	2.5cm Sinus Pack	2/Pack, 10/Box
AM20200-L	2.5cm Sinus Pack, Layered	2/Pack, 10/Box
AM20201-C	2.5cm Sinus Pack Slim	2/Pack, 10/Box
AM20202-C	2.5cm Sinus Pack, Ultra Slim	2/Pack, 10/Box
AM20210-C	3.5cm Sinus Pack	2/Pack, 10/Box
AM20210-L	3.5cm Sinus Pack, Layered	2/Pack, 10/Box
AM20211-C	3.5cm Sinus Pack, Slim	2/Pack, 10/Box
AM20212-C	3.5cm Sinus Pack, Ultra Slim	2/Pack, 10/Box

Epistaxis Products

The Ultracell® brand of epistaxis sponges provides highly absorbent treatment of nosebleeds.

- Fluid absorption causes the packing to expand producing a tamponading effect.
- Available in a variety of sizes: 3.0cm anterior to 10.0cm posterior
- All are strung for safety
- Ideal for use in emergency rooms, clinics, and doctor's offices
- Latex-free
- Single-use



Part No.	Description	Quantity
AM10140-C	10.0cm Epistaxis Pack	1/Pack, 10/Box
AM10140-L	10.0cm Epistaxis Pack, Layered	1/Pack, 10/Box
AM10141-C	10.0cm Epistaxis Slim Pack	1/Pack, 10/Box
AM10145-C	8.0cm Epistaxis Pack	1/Pack, 10/Box
AM10146-C	8.0cm Epistaxis Slim Pack	1/Pack, 10/Box
AM10150-C	5.5cm Epistaxis Pack	1/Pack, 10/Box
AM10150-L	5.5cm Epistaxis Pack, Layered	1/Pack, 10/Box
AM10151-C	5.5cm Epistaxis Slim Pack	1/Pack, 10/Box
AM10160-C	3.0cm Anterior Epistaxis Pack	1/Pack, 10/Box

Otology Products

- American Surgical provides both ear wicks and ear packs.
- Much smaller, compressed form allows easy insertion into swollen ear canal with minimal patient discomfort
 - Pediatric sizes available



Ear Wicks

- Perfect for use with treatments for otitis externa
- Help maintain antibiotics in the ear canal
- Also available fenestrated
- Latex-free
- Single-use

Part No.	Description	Quantity
AM 30300-C	Ear Wick (9mm x 15mm)	1/Pack, 50/Box
AM 30301-C	Pediatric Ear Wick (7mm x 12mm)	1/Pack, 50/Box
AM 30305-C	Ear Wick, Fenestrated (9mm x 15mm)	1/Pack, 50/Box
AM 30306-C	Pediatric Ear Wick, Fenestrated (7mm x 12mm)	1/Pack, 50/Box
AM 30310-C	Ear Wick, Long (9mm x 24mm)	1/Pack, 10/Box
AM 30311-C	Pediatric Ear Wick, Long (7mm x 20mm)	1/Pack, 10/Box
AM 30315-C	Ear Wick, Fenestrated, Long (9mm x 24mm)	1/Pack, 10/Box
AM 30316-C	Pediatric Ear Wick, Fenestrated, Long (7mm x 20mm)	1/Pack, 10/Box




Ear Packs

- Used following otic surgery
- Provides post-surgical packing and dressing of ear canal
- Latex-free
- Single-use

Part No.	Description	Quantity
AM 30320-C	Ear Pack (9mm x 15mm)	1/Pack, 10/Box
AM 30330-C	Ear Pack Large (12mm x 15mm)	1/Pack, 10/Box
AM 30340-C	Ear Pack Large, Long (12mm x 24mm) No string	1/Pack, 10/Box

Otolaryngology



Micro Spheres

For tissue manipulation during ear surgery. Latex-free and intended for single-use.

Part No.	Description	Quantity
AM 30380-C	3mm Micro Sphere	10/Pack, 20/Box
AM 30381-C	4mm Micro Sphere	10/Pack, 20/Box
AM 30382-C	5mm Micro Sphere	10/Pack, 20/Box

Specialty Products

American Surgical also provides specialty PVA products. Included are PVA wipes used in the care and cleaning of surgical instruments. We also distribute a PVA surgical headband. Wipe and cleaning block intended for single-use. Headband can be washed and re-used. All are latex-free.



Part No.	Description	Quantity
AM 40462	Diamond Knife Cleaning Block	1/Pack, 20/Box
AM 90900	Instrument Wipe (2mm)(3"x3"x2mm)	1/Pack, 20/Box
AM 90910	Instrument Wipe (1mm)	1/Pack, 20/Box
AM 90950	Headband	1/Pack, 20/Box

ULTRA POINT

Ophthalmic Knives & Blades

ANGLED KNIVES AND BLADES



Full Handle Knife	Screw-in Blade	Description	Package Contents
STAB AND PARACENTESIS			
AM961501	AM961502	15°	6/Box
AM922201	AM922202	22.5°	6/Box
AM923001	AM923002	30°	6/Box
AM924501	AM924502	45°	6/Box

Full Handle Knife	Screw-in Blade	Description	Color Code	Package Contents
CONTROLLED DEPTH				
AM7512	AM7511	15°, 1.5mm Depth	Orange	6/Box
AM7514	AM7513	15°, 3.0mm Depth	Lt. Blue	6/Box
AM7516	AM7515	15°, 5.0mm Depth	Green	6/Box
AM7518**	AM7517**	15°, 7.0mm Depth	Black	6/Box
AM7523	AM7522	22.5°, 4.0mm Depth	Grey	6/Box
AM7531	AM7530	30°, 3.5mm Depth	Black	6/Box
AM7536**	AM7535**	30°, 5.0mm Depth	Black	6/Box
AM7546**	AM7545**	45°, 3.0mm Depth	Tan	6/Box
AM7551**	AM7550**	45°, 5.0mm Depth	Tan	6/Box

Handles are color coded for easy identification of size and depth.

**Minimum order required. Please allow 4 weeks for delivery.

CRESCENT KNIVES AND BLADES



Full Handle Knife	Screw-in Blade	Description	Package Contents
CRESCENT ANGLED - Round Stock			
AM960022	AM960012	Angled, Slim Design	6/Box


Full Handle Knife	Screw-in Blade	Description	Package Contents
CRESCENT STRAIGHT - Round Stock			
AM960021	AM960011	Straight, Slim Design	6/Box

SLIT KNIVES AND BLADES



Single Bevel Knife	Double Bevel Knife	Description	Package Contents
PHACO SLIT ANGLED - Round Stock			
AM962261		2.2mm Angled	6/Box
AM962361		2.3mm Angled	6/Box
AM962561	AM962568	2.5mm Angled	6/Box
AM962661	AM962668	2.65mm Angled	6/Box
AM962751		2.75mm Angled	6/Box
AM962861	AM962868	2.8mm Angled	6/Box
AM962851		2.85mm Angled	6/Box
AM962961	AM962968	3.0mm Angled	6/Box
AM963261	AM963268	3.2mm Angled	6/Box
AM963561		3.5mm Angled	6/Box
AM965061		5.0mm Angled	6/Box

MVR KNIVES



Full Handle Knife	Description	Package Contents
MVR STRAIGHT		
AM965561	19G (1.1mm), Straight	6/Box
AM965560	20G (.90mm), Straight	6/Box
AM965562	23G (.60mm), Straight	6/Box
AM965564	25G (.50mm), Straight	6/Box









Full Handle Knife	Description	Package Contents
MVR ANGLED		
AM965566	19G (1.1mm), Angled	6/Box
AM965565	20G (.90mm), Angled	6/Box
AM965567	23G (.60mm), Angled	6/Box
AM965569	25G (.50mm), Angled	6/Box

IMPLANT KNIVES

Full Handle Knife	Description	Package Contents
IMPLANT - Rounded Stock		
AM903064	3.0mm Angled, Rounded Tip	6/Box
AM903264	3.2mm Angled, Rounded Tip	6/Box
AM903564	3.5mm Angled, Rounded Tip	6/Box
AM903864	3.8mm Angled, Rounded Tip	6/Box
AM904061	4.0mm Angled, Rounded Tip	6/Box
AM904161	4.1mm Angled, Rounded Tip	6/Box
AM905161	5.1mm Angled, Rounded Tip	6/Box
AM905261**	5.2mm Angled, Rounded Tip	6/Box

**Minimum order required. Please allow 4 weeks for delivery.

MINIATURE BLADES AND KNIVES

			
AM6100 Chisel Small Shaft 12/box	AM6200 Chisel 12/box	AM6300 Spear Tip 6/box	AM6400 Round Tip, Sharp One Side AM4764 Knife Version 12/box 6/box
			
AM6500 Lance Tip 12/box	AM6700 Gradual Curved Tip box	AM6900 Round Tip, Sharp All Around AM4769 Knife Version 12/ 6/box 6/box	AM6910 Round Tip, Sharp All Around, Angled 6/box

SCLERAL KNIVES AND BLADES



Full Handle Knife	Screw-in Blade	Description	Package Contents
SCLERAL			
AM5710	AM5700	Standard USA Size	6/Box
AM5710S	AM5700S	Small European Style	6/Box

OTHER KNIVES AND BLADES



Full Handle Knife	Screw-in Blade	Description	Package Contents
PRECISE DEPTH INCISION			
AM5485	AM5435	.35mm Depth	6/Box
AM5486	AM5440	.40mm Depth	6/Box
AM5490	AM5450	.50mm Depth	6/Box
AM5491	AM5460	.60mm Depth	6/Box

Screw-in Blade	Description	Package Contents
LAMELLAR		
AM6600	Angled, Bevel Up	6/Box
AM6610	Straight	6/Box

**Minimum order required. Please allow 4 weeks for delivery.








CORNEAL PRODUCTS - TREPHINES



Code	Description	Package Contents
TRADITIONAL STYLE		
AM330500	5.00mm	1/Box
AM330550	5.50mm	1/Box
AM330600	6.00mm	1/Box
AM330625	6.25mm	1/Box
AM330650	6.50mm	1/Box
AM330675	6.75mm	1/Box
AM330700	7.00mm	1/Box
AM330725	7.25mm	1/Box
AM330750	7.50mm	1/Box
AM330775	7.75mm	1/Box
AM330785	7.85mm	1/Box
AM330800	8.00mm	1/Box
AM330825	8.25mm	1/Box
AM330850	8.50mm	1/Box
AM330875	8.75mm	1/Box
AM330900	9.00mm	1/Box
AM331000	10.00mm	1/Box
AM331100	11.00mm	1/Box
AM331200	12.00mm	1/Box
AM331300	13.00mm	1/Box
EYE BANK STYLE		
AM331700**	17.00mm	1/Box
AM331800**	18.00mm	1/Box
AM331900**	19.00mm	1/Box

Other sizes available upon request. Minimum order required. Please allow 6-8 weeks for delivery.

ANESTHESIA



Code	Description	Package Contents
------	-------------	------------------

RETROBULBAR, Atkinson Tip

AM111638	23G x 1 1/2 in. (.60 x 38mm)	10/Box
AM111637	25G x 1 1/2 in. (.50 x 38mm)	10/Box

Code	Description	Package Contents
------	-------------	------------------

RETROBULBAR CURVED, Atkinson Tip

AM111284	25G x 1 1/2 in. (.50 x 38mm)	10/Box
----------	------------------------------	--------

Code	Description	Package Contents
------	-------------	------------------

RETROBULBAR, Sharp Tip

AM111278	25G x 1 in. (.50 x 25mm)	10/Box
AM111275	25G x 1 1/2 in. (.50 x 38mm)	10/Box
AM111282	30G x 1/2 in. (.30 x 12.5mm)	10/Box

Code	Description	Package Contents
------	-------------	------------------

PERIBULBAR, Atkinson Tip




AM115144	23G x 3/4 in. (.60 x 22mm)	10/Box
AM115027	25G x 3/4 in. (.50 x 22mm)	10/Box
AM115145	27G x 3/4 in. (.40 x 22mm)	10/Box
AM111641	23G x 1 1/2 in. (.60 x 32mm)	10/Box
AM111642	25G x 1 1/2 in. (.50 x 32mm)	10/Box
AM115109	27G x 1 1/2 in. (.40 x 32mm)	10/Box

Code	Description	Package Contents
------	-------------	------------------

SUB-TENON

AM125176	19G x 3/4 in. (1.1 x 22mm)	10/Box
----------	----------------------------	--------

ANTERIOR CHAMBER MAINTAINERS



Code	Description	Package Contents
------	-------------	------------------

RIBBED OPEN TIP











AM125061	20G (.90)	3/Box
AM125062	23G (.60)	3/Box

Code	Description	Package Contents
------	-------------	------------------

SMOOTH OPEN TIP

AM125163	20G (.90)	10/Box
AM121628	27G (.40)	10/Box

ANTERIOR CHAMBER (RYCROFT)



Code	Description	Package Contents
ANGLED TIP 4mm		
AM121281	20G (.90mm)	10/Box
AM125025	23G (.60mm)	10/Box
AM121639	25G (.50mm)	10/Box
AM101280	27G (.40mm)	10/Box
AM121273	30G (.30mm)	10/Box


Code	Description	Package Contents
ANGLED TIP 6mm		
AM121287	20G (.90mm)	10/Box
AM121636	25G (.50mm)	10/Box
AM121286	27G (.40mm)	10/Box
AM121276	30G (.30mm)	10/Box

Code	Description	Package Contents
ANGLED TIP 7mm		
AM121291	20G (.90mm)	10/Box
AM121637	25G (.50mm)	10/Box
AM121290	27G (.40mm)	10/Box
AM121277	30G (.30mm)	10/Box

Code	Description	Package Contents
ANGLED TIP 8mm (LASIK)		
AM121296	20G (.90mm)	10/Box
AM121638	25G (.50mm)	10/Box
AM121293	27G (.40mm)	10/Box
AM121278	30G (.30mm)	10/Box










Code	Description	Package Contents
ANGLED TIP 10mm (LASIK)		
AM121280	25G (.50mm)	10/Box
AM121274	27G (.40mm)	10/Box

CYSTOTOMES



Code	Description	Package Contents
FORMED		
AM131599	23G x $\frac{7}{8}$ in. (.60 x 16mm)	10/Box
AM131610	25G x $\frac{7}{8}$ in. (.50 x 16mm)	10/Box
AM131618	27G x $\frac{7}{8}$ in. (.40 x 16mm)	10/Box
AM131630	30G x $\frac{7}{8}$ in. (.30 x 16mm)	10/Box

CYSTOTOMES



Code	Description	Package Contents
STRAIGHT		
AM131594	23G x $\frac{3}{8}$ in. (.60 x 16mm)	10/Box
AM131617	25G x $\frac{3}{8}$ in. (.50 x 16mm)	10/Box
AM135017	27G x $\frac{3}{8}$ in. (.40 x 16mm)	10/Box
AM131631	30G x $\frac{3}{8}$ in. (.30 x 16mm)	10/Box

Code	Description	Package Contents
ANGLED 45°		
AM135049	25G x $\frac{3}{8}$ in. (.50 x 16mm)	10/Box
AM135093	27G x $\frac{3}{8}$ in. (.40 x 16mm)	10/Box

Code	Description	Package Contents
REVERSE BEND		
AM135115	25G x $\frac{3}{8}$ in. (.50 x 16mm)	10/Box

Code	Description	Package Contents
ENDOCAPSULAR (Joyce)		
AM135039	27G x $\frac{3}{8}$ in. (.40 x 16mm)	10/Box

Code	Description	Package Contents
SHORT		
AM135004	25G x $\frac{3}{8}$ in. (.50 x 16mm)	10/Box
AM135008	27G x $\frac{3}{8}$ in. (.40 x 16mm)	10/Box



Code	Description	Package Contents
SIDE-CUT		
AM131602	25G x $\frac{3}{8}$ in. (.50 x 16mm)	10/Box
AM131609	27G x $\frac{3}{8}$ in. (.40 x 16mm)	10/Box

Code	Description	Package Contents
PUSH-PULL, Formed		
AM135113	25G x $\frac{3}{8}$ in. (.50 x 16mm)	10/Box
AM135042	27G x $\frac{3}{8}$ in. (.40 x 16mm)	10/Box


Code	Description	Package Contents
PUSH-PULL, Angled 45°		
AM131611	25G x $\frac{3}{8}$ in. (.50 x 16mm)	10/Box

Code	Description	Package Contents
SHARP TIP		
AM135010	22G x 1 in. (.70 x 25mm)	10/Box



HYDRODISSECTORS




Code	Description	Package Contents
FLAT TIP, Angled 7mm		
AM125024	23G (.60mm)	10/Box
AM125036	25G (.50mm)	10/Box
AM125162	27G (.40mm)	10/Box
AM101592	30G (.30mm)	10/Box




Code	Description	Package Contents
FLAT TIP, Angled 8mm		
AM125023	23G (.60mm)	10/Box
AM125037	25G (.50mm)	10/Box
AM125158	27G (.40mm)	10/Box
AM121592	30G (.30mm)	10/Box





Code	Description	Package Contents
HORIZONTAL FLAT TIP, Curved		
AM125155	27G (.40mm)	10/Box





Code	Description	Package Contents
VERTICAL FLAT TIP, Curved		
AM125099	27G (.40mm)	10/Box




Code	Description	Package Contents
CORTICAL CLEAVING, Flat Tip, Angled 11mm		
AM105158	27G (.40mm)	10/Box



Code	Description	Package Contents
CORTICAL CLEAVING, Round Tip, Angled 11mm		
AM125157	27G (.40mm)	10/Box



Code	Description	Package Contents
CORTICAL CLEAVING, Round Tip, Angled 12mm		
AM125197	30G (.30mm)	10/Box

















Code	Description	Package Contents
J-SHAPED, Vertical Flat Tip		
AM125102	25G (.50mm)	10/Box
AM125029	27G (.40mm)	10/Box

FLAT STOCK KNIVES AND BLADES

Screw-in Blade Full Handle Knife Description Pkg Contents


MYRINGOTOMY

Blades require Flat Stock Handle #AM1310 (see below).

AM7100	AM7100H	Lance - Adult Size	6/Box				AM7100H
AM7110		Lance Reverse - Adult Size	6/Box				
AM7121	AM7124	Lance - Juvenile Size	6/Box				AM7124
AM7120	AM7122	Spear - Juvenile Size	6/Box				AM7122
AM7300	AM7300H	Sickle	6/Box				AM7300H

TYMPANOPLASTY


AM4451		Angled 50°, 2mm wide	6/Box
AM7200		Angled 60°, 2.5mm wide	6/Box
AM7210		Straight, 2.5mm wide	6/Box



NEURO


Blade lengths approximately 13.5cm.

AM0110		Sickle, Offset 15°	3/Box
AM0410		Sickle, Bayonet	3/Box
AM0120		Arachnoid, Offset 15°	3/Box
AM0420		Arachnoid, Bayonet	3/Box



SEPTOTOME

AM6001		Angled 90°	6/Box
AM6002		Angled 60°	6/Box



FLAT STOCK HANDLE

Handle accommodates all Flat Stock Blades pictured above.

AM1310	10cm Length with Collet	1/Box
AM1305	Replacement Collet Only	1/Box



ROUND STOCK KNIVES

Full Handle Knife	Description	Pkg Contents
NASAL SINUS SERIES		
AM8900	Spear Tip Angled	4/Box
AM8910	Sickle Tip Curved Left	4/Box
AM8920	Sickle Tip Straight	4/Box
AM8930	Sickle Tip Curved Right	4/Box
AM9610	Uncinate	4/Box



AM8900



AM8910




AM8920



AM8930



AM9610



ROUND STOCK KNIVES AND BLADES

Screw-in Blade	Full Handle Knife	Description	Pkg Contents
MYRINGOTOMY			

Blades require Round Stock Handle #AM1318 (see below).

AM8101	AM8111	Lance Straight - Adult Size	6/Box
AM8102	AM8112	Lance Angled - Adult Size	6/Box



AM8101


AM8102



AM8111

AM8112

AM8121	AM8124	Lance Straight- Juvenile Size	6/Box
AM8123	AM8125	Lance Angled - Juvenile Size	6/Box



AM8121

AM8123



AM8124

AM8125

AM8120	AM8122	Spear Straight	6/Box
AM8119	AM8118	Spear Angled	6/Box



AM8120

AM8119



AM8122

AM8118

AM8130	AM8132	Spear Micro-Point Straight	6/Box
AM8131	AM8133	Spear Micro-Point Angled	6/Box



AM8130

AM8131



AM8132

AM8133

AM8173	Sickle Straight	6/Box
--------	-----------------	-------



AM8173

AM8140	AM8142	Straight Edge - Straight	6/Box
AM8141	AM8143	Straight Edge - Angled	6/Box



AM8140

AM8141



AM8142

AM8143

ROUND STOCK HANDLE

Handle accommodates all Round Stock Blades pictured above.
AM1318 10cm Length

1/Box



AM1318

ROUND STOCK KNIVES AND BLADES

Full Handle Knife Description Pkg Contents

MYRINGOTOMY

AM8324 Bayonet Spear 6/Box



AM8324



Hand-held Blade Description Pkg Contents

EXTENDED REACH MYRINGOTOMY

Longer shaft and finger grip eliminates need for handle. All blade lengths 166mm.

AM8520 Spear Blade 6/Box
(same blade design as AM8120)

AM8521 Lance Blade - Juvenile Size 6/Box
(same blade design as AM8121)

AM8522 Lance Blade - Adult Size 6/Box
(same blade design as AM8101)



AM8520



AM8521



AM8522



Screw-in Blade Description Pkg Contents

TAPERED SERIES

Blades require Tapered Round Stock Handle #AM1382 (see below).

AM8230 Straight Edge - Angled Upward Cut 6/Box

AM8231 Straight Edge - Angled Downward Cut 6/Box

AM8232 Straight Edge - Straight Cut 6/Box

AM8233 Spear Tip Angled 6/Box

AM8234 Spear Tip Straight 6/Box



AM8230



AM8231



AM8232



AM8233



AM8234



TAPERED ROUND STOCK HANDLE

Handle accommodates all Tapered Round Stock Blades pictured

AM1382 10cm Length 1/Box



AM1382

CARDIOVASCULAR

AM4410 DeBakey Blade, 5mm 24/Box

AM4710 DeBakey Blade, 8mm 24/Box

AM4810 DeBakey Blade, 10mm 24/Box



AM4410



AM4710



AM4810

VISCOT Surgical Skin Markers That Last!

- Scrub Resistant Inks
- Durable Tips
- Non-smearing
- Available in sterile & non-sterile



1451-200, 1451-1000 (non-Sterile)

1451SR-100 (Sterile)

Pre-Surgical Mini Marker, Fine Tip. In plastic containers of 100 markers.



Twin Surgical Tip
Marker #1425

1425SR-100 (Sterile) **1425-100** (non-Sterile)

Twin Surgical Skin Marker with Regular Tip one side/Ultra Fine on the opposite side. In plastic containers of 100 markers.



#1436

1436SR-100 (Sterile) **1436-500** (non-Sterile)

Value Surgical Skin Marker. Ultra Fine Tip with Ruler. In plastic containers of 100 markers.



1437SR-100 (Sterile) **1437-500** (non-Sterile)

Value Surgical Skin Marker Reg Tip with Ruler. In plastic containers of 100 markers.



1450XL-200, 1450XL-1000 (non-Sterile)

MINI XL Pre-Surgical Skin Marker, Fine Tip Non-Sterile. In plastic containers of 100 markers.

Drapes

Durable Disposable Products for X-Ray and Imaging Applications



C-Arm

Item Number	Item Description	Unit Packing	Compatible With
AM1145	9" I.I. / X-Ray Tube Cover <i>APPLICATION: 30"X32" Contoured Bandbag. Sterile, individually pouched.</i>	Box of 60	All 9" Image Intensifiers, All X-Ray Tubes.
AM2145	12" I.I. / X-Ray Tube Cover <i>APPLICATION: 36"X36" Contoured Bandbag. Sterile, individually pouched.</i>	Box of 30	All 12" Image Intensifiers, All X-Ray Tubes.
AM4145	9" I.I. Universal Halfbag Cover <i>APPLICATION: 36"X80" Designed to fit the I.I. and midway down the C-Arm. Sterile, individually pouched.</i>	Box of 20	All C-Arms with 9" Image Intensifier.
AM7245	Super C-Arm Drape for OEC 9600/9800 <i>APPLICATION: Consists of C-Arm Drape with 5 mounting clips. Sterile, individually pouched.</i>	Box of 30	OEC 9600/9800, Siemens ARCADIS Avantic, Siemens SIREMOBIL Compact L, Siemens ARCADIS Varic, Philips BV Series, Ziehm.
AM9245	12" I.I. C-Arm Kit with Clips, Fits OEC 9600/9800" <i>APPLICATION: Consists of (AM7245)C-Arm Drape with 5 mounting clips, (AMT756)I.I. Cover, (AMT756)X-Ray Tube Cover and a (AM1245)Footswitch Cover, Sterile, individually pouched.</i>	Box of 20	OEC 9600/9800, Siemens ARCADIS Avantic, Siemens SIREMOBIL Compact L, Siemens ARCADIS Varic, Philips BV Series, Ziehm.
AM4345	Universal Mobile C-arm Drape 42" x 74" <i>APPLICATION: Clear poly, 42"x74" with poly straps. Sterile, individually pouched.</i>	Box of 20	All C-Arms with 9" or 12" Image Intensifier.
AM0645	Siemens Arcadis Orbic 3D - Complete C-Arm Kit <i>APPLICATION: Kit consists of C-Arm/I.I. Cover X-Ray Tube Cover, Footswitch Cover and Keyboard Cover. Sterile, individually pouched.</i>	Box of 20	Siemens SIREMOBIL ISO-C & Siemens ARCADIS Orbic.
AM1645	Siemens Arcadis Orbic 3D/C-Arm Drape Only <i>APPLICATION: Consists of C-Arm/I.I. Cover. Sterile, individually pouched.</i>	Box of 20	Siemens SIREMOBIL ISO-C & Siemens ARCADIS Orbic
AM4645	Siemens Arcadis Orbic 3D/C-Arm Drape & X-Ray Tube Cover <i>APPLICATION: Consists of (AM1645)C-Arm/I.I. Cover and X-Ray Tube Cover. Sterile, individually pouched.</i>	Box of 20	Siemens SIREMOBIL ISO-C & Siemens ARCADIS Orbic.

Drapes

Durable Disposable Products for Imaging Accesories



Mini C-Arm

Item Number	Item Description	Unit Packing	Compatible With
AM2045	Mini C-Arm Kit <i>APPLICATION: Kit consists of a (AM2045) Mini View C-Arm Drape, (AM5045) Keyboard Cover and a (AM1245) Footswitch cover. Sterile, individually pouched.</i>	Box of 20	Fits OEC Mini 6600/6800
AM3045	Universal Mini C-Arm Drape <i>APPLICATION: Single piece drape to cover I.I. X-Ray Tube and the C-Arm. Sterile, individually pouched.</i>	Box of 25	Fits All Mini C-Arms

Dome & Band Bags

Item Number	Item Description	Unit Packing	Application
AM756	36" X 36" Band Bag	Box of 100	Generic Equipment Cover
AM766	30" X 36" Band Bag	Box of 100	Generic Equipment Cover
AM657	54" expanded / 27" depth Dome Bag	Box of 50	Generic Equipment Cover
AM128	44" expanded / 22" depth Dome Bag	Box of 50	Generic Equipment Cover
AM728	40" X 40" Band Bag	Box of 100	Generic Equipment Cover



Footswitch Covers



Item Number	Item Description	Unit Packing	Compatible With
AM6045	Footswitch Cover <i>APPLICATION: 17" x 15" OEC Gray Footswitch. Fits with cinchlock</i>	Box of 100	Fits OEC Gray Footswitch
AM7045	Footswitch Cover <i>APPLICATION: 10" x 8" OEC Black Footswitch. Fits with cinchlock</i>	Box of 100	Fits OEC Black Footswitch
AM1245	Wide Footswitch Cover <i>APPLICATION: 24" x 20" OEC UroView 2600 Wide Footswitch. Fits with cinchlock</i>	Box of 100	Fits OEC UroView 2600 Wide Footswitch
AM2245	Wide Footswitch Cover <i>APPLICATION: 20" x 36" OEC UroView 2800 Wide Footswitch. Fits with cinchlock</i>	Box of 100	Fits OEC UroView 2800 Wide Footswitch




Drapes

Durable Disposable Products for Ultrasound and Urology



Ultrasound Covers

Item Number	Item Description	Unit Packing	Compatible With	Application
AM7845	Ultrasound Cover Kit	Box of 20	Ultrasound Cover Kit 8.9cm x 61cm	Sterile
AM8845	Ultrasound Cover Kit	Box of 20	Ultrasound Cover Kit 14cm x 61cm	Sterile
AM9845	Ultrasound Cover Kit	Box of 20	Ultrasound Cover Kit 17.8 cm x 61cm	Sterile
AM0945	Ultrasound Cover Kit	Box of 20	Ultrasound Cover Kit 22.9 cm x 61cm	Sterile
AM1945	Ultrasound Cover Kit	Box of 20	Ultrasound Cover Kit Tapered 12.7-5cm x 61cm	Sterile
AM2945	Ultrasound Cover Kit	Box of 20	Ultrasound Cover Kit 14cm x 91cm	Sterile
AM3945	Ultrasound Cover Kit	Box of 20	Ultrasound Cover Kit 15.2cm x 122cm	Sterile
AM5945	Ultrasound Cover Kit	Box of 20	Ultrasound Cover Kit Tapered 12.7cm-5cm x 244cm	Sterile
AM7945	Ultrasound Cover Kit	Box of 100	Ultrasound Cover Kit 8.9 cm x 61cm	Non - Sterile



Urology Covers

Item Number	Item Description	Unit Packing	Compatible With
AM0005	Square Drain Bag w/ Metal Tabs <i>APPLICATION: All Square Urology Frames fit the following products: GE/OEC Uroview 2000 and GE/ OEC Uroview 2500, Liebel Flarshheim/ Mallinckrodt Hydrovision ED, DR, and HP Urology Systems with the collapsible bracket, fixed bracket, or square bracket.</i>	Box of 20	Fits All Square Frames
AM0105	Urology Drain Collection Bag <i>APPLICATION: All Urology Drain Bags. Collects waste fluid when floor drain or suction is unavailable. Internal baffle system enables container to fill evenly as powder converts liquid to gel. Built in handles to transport unit for disposal. 5 gallon capacity.</i>	Box of 10	
AM3135	Uro-Catcher Bag <i>APPLICATION: Fits all hoop frames. Under buttox flap directs fluids into the drain bag, tissue filter, and flexible hose that extends 6 feet. Sterile, individually pouched.</i>	Box of 20	Fits All Hoop Style Frames
AM6145	Urology Disposables Pack <i>APPLICATION: Includes (AM9145) Drain Bag, (AM0245) Hand Control Cover, (AM1245) 24"x20" Footswitch Cover. Sterile, individually pouched.</i>	Box of 20	Fits OEC UroView 2600/2800
AM7145	Square Drain Bag w/ Plastic Tabs <i>APPLICATION: All Square Urology Frames fit the following products: GE/OEC Uroview 2000 and GE/ OEC Uroview 2500, Liebel Flarshheim/ Mallinckrodt Hydrovision ED, DR, and HP Urology Systems with the collapsible bracket, fixed bracket, or square bracket.</i>	Box of 20	Fits All Square Frames
AM9145	Hoop Drain Bag <i>APPLICATION: Fits GE/OEC Uroview 2600 and GE/OEC Uroview 2800 tables. Heavy gauge vinyl strip holds bag snugly to the foot of the table, Mesh screen, flexible hose extends 6 feet. Sterile, individually packaged.</i>	Box of 20	Fits OEC UroView 2600/2800
AM2445	Urology Drain Bag for Siemens Urooskop System <i>APPLICATION: Siemens Urooskop & Dornier Systems</i>	Box of 20	
AM0845	BAG, URO-CATCHER POLY <i>APPLICATION: All Hoop Style Frames. Hoop Style Frames are flexible Rubber/Metal frames that clamp on to the rail of an operating room table</i>	Box of 20	

Drapes

Durable Disposable Products for Microscopes, Keyboards and Controllers



Microscope Covers

Item Number	Item Description	Unit Packing	Compatible With
AM562144	Ziess MD Microscope Drape <i>APPLICATION: Strong drapable polyethylene film, monocular or stereo observer tubes, including straps to gather excess poly and re-attachable fiber-free ocular closures.</i>	Box of 10	44" wide x 120.0" long fits 65 mm lens
AM562124	Ziess OPMI Microscope Drape	Box of 10	42.0" wide x 120.0" long fits 65mm lens
AM050644	Ziess OPMI Microscope Drape <i>APPLICATION: Contoured to reduce unnecessary excess drape</i>	Box of 10	44.0" wide x 60.0" long fits 50mm lens
AM865145	Leica Wild Microscope Drape	Box of 10	54.0" wide x 150.0" long fits 68mm lens
AM1345	Video Camera/Laser Arm Drape with Perforated Tip <i>APPLICATION: Camera drapes are telescope folded with cuff to allow for easy insertion. All have double-sided tape for gathering material and secure closure. Latex free with easy to follow instruction labels. Sterile, individually pouched.</i>	Box of 25	Camera Drape, 10" x 80" clear polyethylene with perforated tip
AM2345	Camera Drape for use with Video System	Box of 25	Camera Drape, 5" x 96" clear polyethylene with perforated tip
AM3345	Universal Video Camera/Laser Drape	Box of 25	Blue Tip Camera Drape, 7" x 96" clear polyethylene with elasticized fenestration to provide a secure fit around the hand piece
AM0051	Camera Cover with Blue PE <i>APPLICATION: Blue poly, Sterile, individually pouched</i>	Box of 20	Blue Camera Cover 5" x 96"
AM1051	Blue Probe Cover <i>APPLICATION: Blue poly, Sterile, individually pouched</i>	Box of 20	Blue Ultrasound Probe Cover 5" x 72"
AM4051	Camera Cover Kit <i>APPLICATION: Clear poly, Sterile, individually pouched</i>	Box of 20	Clear Camera Cover 7" x 96"

Miscellaneous Covers



Item Number	Item Description	Unit Packing	Application
AM5405	Mini C Keyboard Cover	Box of 25	Most O.R. Keyboards
AM3145	IR Hand Control Cover, 7" x 13"	Box of 25	7"x13" clear cover. Sterile, individually pouched.
AM207	44" X 70" Set-Up Cover	Box of 25	General Set Up Cover. With Caution Label. Sterile, Individually packaged.
AM407	Slush Drape, 54" X 54"	Box of 25	General Set Up Cover. Sterile, Individually packaged.
AM807	Mayo Stand Cover 24" X 58"	Box of 25	General Mayo Stand Cover.
AM939	Set-Up Cover 36" X 84"	Box of 100	Sterile, individually packaged. General Mayo Stand Cover

Shielding Products

Durable, Flexible Protective Aprons and Skirts
Additional Styles & Colors Available

Your American Surgical specialist will assist you in ordering a stock apron or customizing your very own model. We offer a variety of styles, patterns, and colors including pocket designs, embroidery options, and custom binding. Our most popular styles are listed below with regular lead, TL lite ply, or non-lead. Contact us for a quote and to hear the many options we offer for your protection needs.

Flex Back Aprons




Model - FB Application

Flex Back Apron features elastic back panels for added comfort and support hook and loop fastening flaps that wrap around waist for an adjustable fit. Shoulder pads included.


Two Piece Flex Back Vest Skirts

Model - FBVS Application



Flex Back Vest features wrap around protection with flex back support, adjustable hook and flaps at the waist and shoulder for added comfort and support, especially during long procedures. Shoulder Pads included.

Standard Coat Aprons




Model - S Application

Standard Coat Apron features open back with choice of durable 1" tie straps, quick release belt, or hook and loop fastener that crosses the back and secures in front for a snug and comfortable fit. Shoulder pads included. Protection - .5mm PB Equivalency

Adjust-A-Fit Aprons

Model - AF Application



Adjust - A - Fit Apron features two hook and loop fastening flaps at the waist, which aid in relieving weight off the shoulders, and includes adjustable padded shoulders. Its comfortable fit is attained by varying the position of the flaps. Shoulder pads included.

Surgical Drop Away




Model - SDA Application

Surgical Drop Away/Flex Back Apron offers top of the shoulder hook and loop fastening flaps for an adjustable fit. Features elastic back panels for added comfort and support. For quick removal, apron releases at shoulder or waist. Shoulder pads included.

Surgical Drop Off Aprons

Model - SDO Application



Surgical Drop Off Aprons. Features back criss-cross flaps with front hook and loop closure that assure easy removal from both over or under the surgical gown without disturbing the sterile field. Shoulder pads included.

SHIELDING INTERNATIONAL™ Offers (3) Choices for Protection
Choose from Regular (Lead, optional TL Lite Ply, or environmentally safe NL Non Lead option

Standard Protection Regular Lead	Benefits Most Economical 100 % Recyclable	Protection Value Provides 0.5mm Pb Equivalent Protection @ 100 kVp - 96 % Attenuative
Optional Protection TL Lite Ply	Benefits 20 % Lighter than Regular Lead 100 % Recyclable	Protection Value Provides 0.5mm Pb Equivalent Protection @ 100 kVp - 96 % Attenuative
Optional Protection NL Non Lead	Benefits Lighter than Regular Lead Lead Free & Environmentally Friendly	Protection Value Provides 0.5mm Pb Equivalent Protection @ 90 kVp - 99% Attenuative


Product Conformance: In compliance with ASTM test method F2507-06. Core protection material testing performed by Health Physics Northwest.

2 Year Limited Protective Apparel Warranty!

SHIELDING INTERNATIONAL™ warrants Regular Lead and TL Lite Ply x-ray protective apparel from valid material defects and workmanship for a period of (2) years from date of purchase.
(NL, Non Lead apparel is limited to a 90 day warranty.)



Wrap, Demi & Half Aprons



Application

The three quarter wrap apron is especially designed for the Fluoro Field. It features 3/4 coverage for added protection and increased mobility. Adjustable padded shoulder flaps with elastic back panels for added comfort. 0.5mm Pb Eq standard protection.

Demi Apron - Coated spring ring closure or QR buckle. Available in sizes S-XL. Also available as a set of 4 with rack. 0.5mm Pb Eq Protection. (not pictured)
Half Apron - Closures: Ties, Hook & Loop or QR buckle. Ideal for front or back protection. Available in assorted sizes. 0.5mm Pb Eq Protection. (pictured)



Application

Standard Protection Accesories



STANDARD

AM600S Small
AM600L Large



NEW STYLE

AM611NSS Small
AM611NSL Large



ULTRA

AM621S Small
AM621L Large



BUCKLE

AM631QRS Small
AM631QRL Large



COMFORT FIT

AM641CFS Small
AM641CFL Large



Gonadal & Ovarian Belt

AM822BBV 4Shields

Protective Eyewear



Ultra Light Viper - AMULV99 & AMULV99-S (small)

Contours around face for added protection. Front lenses 0.75 mm Pb Eq optical lead glass. Optional side shields. Available in prescription. Colors: Silver, Blue & Red

Sport Flex - AMSP92 & AMSF92-S (small)

Enhanced wrap around fit. Front lenses 0.75 mm PB Eq optical lead glass. Optional side shields. Available in prescription. Colors: Black, Silver & White



American Surgical

Specialties Company



www.americansurg.com

Copyright 2010
American Surgical Specialties Company

